


Temporary Snap Abutment, Engaging/Non-Engaging

Instructions for use



Important: Please read.

Disclaimer of liability:

This product is part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Nobel Biocare. Non-recommended use of products made by third parties in conjunction with Nobel Biocare products will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare products has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The user is also obliged to study the latest developments in regard to this Nobel Biocare product and its applications regularly. In cases of doubt, the user has to contact Nobel Biocare. Since the utilization of this product is under the control of the user, they are his/her responsibility. Nobel Biocare does not assume any liability whatsoever for damage arising thereof. Please note that some products detailed in this Instructions for Use may not be regulatory cleared, released or licensed for sale in all markets.

Description:

A premanufactured dental implant abutment directly connected to the endosseous dental implant intended for use as a temporary aid in prosthetic rehabilitation.

Temporary Snap Abutment Engaging and Temporary Snap Abutment Non-Engaging can be used with internal conical connection implants NobelActive®, NobelReplace® CC and NobelParallel™ CC.

Intended use:

Dental implant abutments are intended to be used in the upper or lower jaw and used for supporting tooth replacements to restore chewing function.

Indications:

Temporary Snap Abutment Engaging is indicated for single unit screw-retained temporary restorations.

Temporary Snap Abutment Non-Engaging is indicated for screw-retained multiple temporary restorations, for implants with less than 40° overall divergences to allow path of insertion.

Contraindications:

Temporary Snap Abutments are contraindicated for patients:

- who are medically unfit for an oral surgical procedure.
- in whom adequate sizes, numbers or desirable positions of implants are not reachable to achieve safe support of functional or eventually parafunctional loads.
- allergic or hypersensitive to commercially pure titanium, titanium alloy Ti-6Al-4V (titanium, aluminum, vanadium).

Cautions:

Pre-operative hard tissue or soft tissue deficits may yield a compromised esthetic result or unfavorable implant angulations.

To secure the treatment outcome it is advised to provide comprehensive regular patient follow up after implant treatment and to inform about appropriate oral hygiene.

All instruments and tooling used in surgery must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components. Special attention has to be given to patients who have local or systemic factors that could interfere with the healing process of either bone or soft tissue or the osseointegration process (e.g., cigarette smoking, poor oral hygiene, uncontrolled diabetes, oro-facial radiotherapy, steroid therapy, infections in the neighboring bone).

Special caution is advised in patients who receive bisphosphonate therapy.

In general, implant placement and prosthetic design must accommodate individual patient conditions. In case of bruxism or unfavorable jaw relationships reappraisal of the treatment option may be considered.

Never exceed the recommended maximum prosthetic tightening torque for the abutment screw (see clinical procedure). Overtightening of abutment may lead to a screw fracture.

Because of the small sizes of the devices, care must be taken that they are not swallowed or aspirated by the patient.

It is strongly recommended that clinicians, new as well as experienced implant users, always go through special training before undertaking a new treatment method. Nobel Biocare offers a wide range of courses for various levels of knowledge and experience. For more info please visit www.nobelbiocare.com.

Working the first time with a colleague, experienced with the new device/treatment method, avoids eventual complications. NobelBiocare has a global network of mentors available for this purpose.

Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment.

The use of Temporary Snap Abutments is limited to 180 days.

Temporary Snap Abutments shall be taken out of occlusion and should not be used for full-arch restoration.

Handling instructions:

Modifications of abutments could be performed using copious water irrigation. Extra-oral modification of abutment is recommended. Use a carborundum disk and carbide bur.

Clinical procedure (Chair-side made provisionals):

1. Connect the Temporary Snap Abutment to the implant and modify the abutment if necessary using copious irrigation.

Note: Until the Temporary Snap Abutment is secured with the abutment screw, care should be exercised that the Temporary Snap Abutment does not detach from the Implant (e.g. through pressure from the tongue).

2. Close the screw access hole.
3. Make a temporary restoration using a pre-fabricated mold with suitable temporary crown and bridge material.
4. Drill a hole through the mold, loosen the screw(s) using a Unigrip™ Screwdriver and remove the restoration.
5. Make final adjustments.
6. Connect the temporary restoration using a Unigrip™ Screwdriver.
7. Tighten abutment to **35Ncm** using Unigrip™ Machine Screwdriver and Manual Torque Wrench Prosthetic. It is recommended to verify the final seating using radiographic imaging.

Warning: Over tightening of abutment screw may lead to a screw fracture.

Materials:

Temporary Snap Abutment Engaging and Non-Engaging: Titanium alloy 90% Ti, 6% Al, 4%V.

Cleaning and sterilization:

Temporary Snap Abutment is delivered non-sterile for single use only. Prior to use clean, disinfect and sterilize the product using the recommended parameters.

Warning: Use of non-sterile components may lead to infection of tissues or infectious diseases.

Caution: Temporary Snap Abutment is a single use product and must not be reprocessed. Reprocessing could cause loss of mechanical, chemical and / or biological characteristics. Re-use could cause cross contamination.

Warning: Do not use device if the packaging has been damaged or previously opened.

Clean the device using manual or automated cleaning, disinfect and dry the device following instructions in Cleaning and Sterilization Guidelines available at www.nobelbiocare.com/sterilization.

Inspect and seal the single device in a pouch and steam sterilize, both the gravity cycle (saturated steam) and pre-vacuum (forced air removal) cycle can be applied, using the following parameters:

Alternative UK: Temperature 134°C (273°F), max 137°C (279°F) for 3 minutes (up to 20 minutes). Dry for 10 minutes in chamber.

Magnetic Resonance (MR) safety information:


The Temporary Snap Abutment has not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Temporary Snap Abutment in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Storage and handling:

The device must be stored and transported in dry conditions in the original packaging at room temperature and not exposed to direct sunlight. Incorrect storage and transportation may influence device characteristics leading to failure.

Disposal:

Disposal of the device shall follow local regulations and environmental requirements, taking different contamination levels into account.

 **Manufacturer:** Nobel Biocare AB, Box 5190, 402 26 Västtra Hamngatan 1, 411 17 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com



Canada license exemption: Please note that not all products may have been licensed in accordance with Canadian law.

Prescription device: Rx only

Symbols Glossary:

The following symbols may be present on the device labeling or in information accompanying the device. Refer to the device labeling or accompanying information for the applicable symbols.



Authorized representative in the European Community



Batch code



Catalogue number



Caution



CE marking



Consult instructions for use



Contains hazardous substances



Contains or presence of phthalate



Date



Date of manufacture



Do not re-sterilize



Do not re-use



Do not use if package is damaged



Double sterile barrier system

Rx Only

For prescription use only



Health care centre or doctor



Keep away from sunlight



Keep dry



symbol.glossary.nobelbiocare.com
ifu.nobelbiocare.com

Link to Online Symbols Glossary and IFU Portal



Magnetic resonance conditional



Manufacturer



Medical device



Non-pyrogenic



Non-sterile



Patient identification



Patient information website



Patient number



Serial number



Single sterile barrier system



Single sterile barrier system with protective packaging inside



Single sterile barrier system with protective packaging outside



Sterilized using ethylene oxide



Sterilized using irradiation



Temperature limit



Tooth number



Upper limit of temperature



Sterilized using steam or dry heat



Unique Device Identifier



Use-by date

EN All rights reserved.

Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Product images in this folder are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product.