creos™

allo.gain

allo.protect

QC-605-F-20 V 2.0
77586 ML

ENGLISH

DESCRIPTION

DONATED HUMAN TISSUE. Tissue grafts are retrieved from deceased human donors. All tissue is retrieved, processed, stored and distributed for use in accordance with the standards of the American Association of Tissue Banks (AATB).

INTENDED USE

Creos allo.gain consists in bone particulates intended for filling and augmenting intraoral and maxillo-facial osseous defects. Creos allo.protect is a biodegradable collagen barrier membrane intended to protect the graft site during healing

DONOR SELECTION

Donor has been determined to be eligible by a Community Tissue Services (CTS) Medical Director at 349 S. Main St., Dayton, OH 45402 based on the results of screening and testing. Screening includes a review of medical and social history, hospital records, infectious disease screening, autopsy report (if performed), and physical exam. Donors are tested and found negative (acceptable) for anti-HIV 1/2, HBsAg, HBV NAT, anti-HBc, anti-HCV, HIV NAT, HCV NAT and syphilis. U.S. Food and Drug Administration (FDA) licensed test kits are used when available. Additional tests, including but not limited to HTLV I/II, may have been performed and were found to be acceptable for transplantation. Communicable disease testing has been performed by a laboratory registered with the FDA and certified under CLIA or equivalent requirements.

Tissue has been processed with Bacitracin and/or Polymyxin B and traces may remain. Demineralized tissue has also been processed with HCl, alcohol, sodium phosphate (monobasic and dibasic) and traces may remain. Tissue has been processed using a proprietary method that has been validated to a Sterility Assurance Level (SAL) of 10⁻⁶. Irradiated tissue is Gamma Irradiated with Cobalt 60.

WARNINGS AND PRECAUTIONS

1. Intended for use in one patient, on a single occasion only.
2. Do not use if package integrity

Distributed by

Nobel Biocare USA, LLC
22715 Savi Ranch Parkway
Yorba Linda, CA 92887
Tel.: +1-800-322-5001
info.usa@nobelbiocare.com

Nobel Biocare Canada Inc.
9133 Leslie Street, Suite 100
Richmond Hill, ON, L4B 4N1
Tel.: +1-800-939-9394
info@canada@nobelbiocare.com
Canadian reg. no.: 100229

Manufactured by/ source establishment

Community Tissue Services
2900 College Drive
Kettering, OH 45420, USA
Tel.: +1-800-684-7783
Canadian reg. no.: 100076

ty has been compromised. Once the user breaks the container seal, the tissue grafts must be transplanted or discarded.

3. Tissue may not be sterilized or re-sterilized.

4. This tissue is intended for use by qualified healthcare specialists such as physicians, dentists, or podiatrists.

5. Although this tissue has been tested and screened for human pathogens, and processed under aseptic conditions, human derived tissue may still transmit infectious agents.

6. Adverse outcomes potentially attributable to this tissue must be reported promptly to Nobel Biocare.

STORAGE CONDITIONS

Freeze-Dried tissue must be stored at ambient temperature or colder.

Tissue may not be stored at liquid nitrogen (LN₂) vapor phase or LN₂ liquid temperatures. It is the responsibility of the Tissue Dispensing Service, Tissue Distribution Intermediary, and/or End User clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to further distribution or transplant.

TISSUE TRACKING

Complete the enclosed Allograft Tracking Form and send it to Nobel Biocare. Federal Regulations (21 CFR 1271.290(b)) and Joint Commission Standards (TS.03.02.01, EP 7) require proper tracking of this tissue. It is the responsibility of the end-user to provide this information, which enables Nobel Biocare to maintain records for the purpose of tracing the tissue post-transplant.

TISSUE PREPARATION

1. Inspect for package integrity and expiration date prior to opening.
2. Tissue contained within double packaging: peel outer package down and aseptically deliver inner package to the sterile field or sterile team member.
3. Tissue in vacuum sealed jars: peel off metal cap and wipe rubber stopper with alcohol or betadine. Using a syringe, inject sufficient saline or air to release vacuum. If vacuum is present, plunger will be drawn down. **DO NOT USE IF VACUUM IS NOT PRESENT.** Remove rubber stop-

per with aid of sterile forceps.

4. Remove tissue from inner package and place in sterile basin and cover with normal saline or isotonic solution of choice. Antibiotics of choice may be added.

5. **IMPORTANT!** Bone Particulate and Pericardium should be reconstituted for 5 minutes. Final determination of allograft reconstitution should be made by the physician prior to use.

6. Tissue should be used as soon as possible after reconstitution. If tissue is to be stored for longer than two hours after reconstitution, it should be refrigerated at 1 to 10°C in an aseptic container for no longer than 24 hours.

7. **IMPORTANT!** Peel away and remove all internal packaging materials from the graft (i.e. gauze or mesh) prior to implantation.

Nobel Biocare and Community Tissue Services make no claims concerning the biological or biomechanical properties of the provided tissue. Nobel Biocare and Community Tissue Services disclaim all liability and responsibility for any misuse of tissue provided for clinical application.

Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of the product. The user is also obliged to study the latest developments in regard to this product and its applications regularly. Since the processing and surgical application of this product is under the control of the user, they are his/her responsibility. Nobel Biocare does not assume any liability whatsoever for damage arising thereof.

Community Tissue Services is accredited by the American Association of Tissue Banks. Community Tissue Services - Center for Tissue Innovation and Research is ISO 13485:2003 certified. Canadian registration number: 100076.

Please contact Nobel Biocare should you require further information.

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