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Table des matières

Informations de compatibilité	3
IFU1001. Implants NobelActive® TiUnite® et NobelActive® TiUltra™	3
IFU1002. Implants NobelParallel™ CC et NobelParallel™ CC TiUltra™	3
IFU1007. Brånemark System® MkIII TiUnite®, NobelSpeedy®, forets hélicoïdaux Twist Drill, forets hélicoïdaux à paliers Twist Step Drill, fraises coniques Counterbore, tarauds Screw Taps.....	4
IFU1008. NobelReplace® Select™ TC.....	6
IFU1010. NobelReplace® Conical Connection, NobelReplace® Conical Connection PMC, NobelReplace® Conical Connection TiUltra™	7
IFU1012. Replace Select™ Tapered TiUnite®, Replace Select™ Tapered PMC	7
IFU1016. Vis de couverture	8
IFU1018. Esthetic Abutment Conical Connection/Brånemark System®/NobelReplace®, Esthetic Abutment Nobel Biocare N1™ TCC	8
IFU1023. Piliers et bases Universal	9
IFU1040. Abutment Release Pins CC	9
IFU1043. Abutment Screw Retrieval Instrumentation.....	9
IFU1048. Forets anatomiques Tapered Drill, forets pour os de forte densité Dense Bone Drills, tarauds Screw Taps	10
IFU1050. Implant NobelZygoma™ 45°.....	10
IFU1051. Implants NobelZygoma™ 0°	10
IFU1057. Vis cliniques, vis de pilier, vis de prothèse	10
IFU1067. Plateau PureSet™ Tray	11
IFU1074. Concept On1™	11
IFU1075. Piliers Multi-unit Abutments et capuchons de cicatrisation Healing Caps	13
IFU1085. Tournevis	13
IFU1086. Transferts d'empreinte.....	13
IFU1087. Nobel Biocare N1™ TiUltra™ TCC System	14

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Nom du produit	Vis de couverture	Guide-implant	Vis	Piliers			Transfert d'empreinte
				de cicatrisation	provisoire	permanent	
Brånemark System® Mk III TiU NP	Cover Screw Brånemark System® NP	Implant Driver NP Brånemark System Implant Driver Wrench Adapter NP	Abutment Screw Multi-unit Angled Prosthetic Screw Multi-unit	Healing Abutment NP	Immediate Temporary Abutment NP	Procera Esthetic Abutment NP	Impression Coping Closed Tray NP
					Plastic Coping Immediate NP	Esthetic Abutment NP	Impression Coping Open Tray NP
					Temporary Abutment NP	15° Esthetic Abutment NP	Impression Coping Open Tray Multi-unit NP
					QuickTemp™ Abutment Conical NP	Snappy™ Abutment NP	Impression Coping Closed Tray Multi-unit NP
					Plastic Coping QuickTemp™ Abutment Conical NP	Universal Base NP	Impression Coping Closed Tray Multi-unit NP
					Temporary Abutment Engaging NP	GoldAdapt™ Engaging NP	Impression Coping Bar Closed Tray Multi-unit NP
					Temporary Abutment Non-Engaging NP	GoldAdapt™ Non-Engaging NP	Temporary Coping Multi-unit Titanium NP
					Temporary Abutment Plastic Engaging NP	Multi-Unit Abutment NP	Temporary Coping Multi-unit Plastic NP
					Temporary Abutment Plastic Non-Engaging NP	Multi-unit Abutment NP 17 °	
Brånemark System® Mk III TiU RP	Cover Screw Brånemark System® RP	Implant Driver RP Brånemark System® Implant Driver Wrench Adapter RP		Healing Abutment RP	Immediate Temporary Abutment RP	Procera Esthetic Abutment RP	Impression Coping Closed Tray RP
					Plastic Coping Immediate Temporary Abutment RP	Esthetic Abutment RP	Impression Coping Open Tray RP
					QuickTemp™ Abutment Conical RP	15° Esthetic Abutment RP	Impression Coping Bar Brånemark System® RP
					Plastic Coping QuickTemp™ Abutment Conical RP	Snappy™ Abutment RP	Impression Coping Open Tray Multi-unit RP
					Temporary Abutment Engaging RP	Universal Base Burn-out Coping RP	Impression Coping Closed Tray Multi-unit RP
					Temporary Abutment Non-Engaging RP	GoldAdapt™ Engaging RP	Impression Coping Bar Closed Tray Multi-unit RP
					Temporary Abutment Plastic Engaging RP	GoldAdapt™ Non-Engaging RP	Temporary Coping Multi-unit Titanium RP
					Temporary Abutment Plastic Non-Engaging RP	Multi-Unit Abutment RP	Temporary Coping Multi-unit Plastic RP
						17° Multi-Unit Abutment RP	
						30° Multi-Unit Abutment RP	
	30° Multi-Unit Abutment Non-Engaging RP						
	Locator Abutment RP						
	Ball Abutment Titanium RP						
Brånemark System® Mk III TiU WP	Cover Screw Brånemark System® WP	Implant Driver WP Brånemark System Implant Driver Wrench Adapter WP	Prosthetic Screw Multi-unit	Healing Abutment WP	Immediate Temporary Abutment WP	Procera Esthetic Abutment WP	Impression Coping Closed Tray WP
					Plastic Coping Immediate Temporary Abutment WP	Esthetic Abutment WP	Impression Coping Open Tray WP
					QuickTemp™ Abutment Conical WP	15° Esthetic Abutment WP	Impression Coping Open Tray Multi-unit WP
					Plastic Coping QuickTemp™ Abutment Conical WP	Snappy™ Abutment WP	Impression Coping Closed Tray Multi-unit WP
					Temporary Abutment Engaging WP	Universal Base Burn-out Coping WP	Impression Coping Bar Closed Tray Multi-unit WP
					Temporary Abutment Non-Engaging WP	GoldAdapt™ Engaging WP	Temporary Coping Multi-unit Titanium WP
					Temporary Abutment Plastic Engaging WP	GoldAdapt™ Non-Engaging WP	Temporary Coping Multi-unit Plastic WP
					Temporary Abutment Plastic Non-Engaging WP	Multi-Unit Abutment WP	
						Locator Abutment WP	

Nom du produit	Vis de couverture	Guide-implant	Vis	Piliers			Transfert d'empreinte
				de cicatrisation	provisoire	permanent	
NobelSpeedy® Shorty NP	Cover Screw Brånemark System® NP	Implant Driver NP Brånemark System® Implant Driver Wrench Adapter NP	Abutment Screw Multi-unit Angled Prosthetic Screw Multi-unit	Healing Abutment NP	Immediate Temporary Abutment NP	Procera Esthetic Abutment NP	Impression Coping Closed Tray NP
Plastic Coping Immediate NP					Esthetic Abutment NP	Impression Coping Open Tray NP	
NobelSpeedy® Groovy NP					Temporary Abutment NP	15° Esthetic Abutment NP	Impression Coping Open Tray Multi-unit NP
					QuickTemp™ Abutment Conical NP	Snappy™ Abutment NP	Impression Coping Closed Tray Multi-unit NP
					Plastic Coping QuickTemp™ Abutment Conical NP	Universal Base NP	Impression Coping Closed Tray Multi-unit NP
					Temporary Abutment Engaging NP	Burn-out Coping NP	Impression Coping Bar Closed Tray Multi-unit NP
					GoldAdapt™ Engaging NP	GoldAdapt™ Engaging NP	Impression Coping Bar Closed Tray Multi-unit NP
					Temporary Abutment Non-Engaging NP	GoldAdapt™ Non-Engaging NP	Temporary Coping Multi-unit Titanium NP
					Temporary Abutment Plastic Engaging NP	Multi-Unit Abutment NP	Temporary Coping Multi-unit Plastic NP
					Temporary Abutment Plastic Non-Engaging NP	Multi-unit Abutment NP 17 °	
						Locator Abutment NP	
						Gold Coping Bar Multi-unit (vis de prothèse incluse) NP	
NobelSpeedy® Shorty RP	Cover Screw Brånemark System® RP	Implant Driver RP Brånemark System Implant Driver Wrench Adapter RP		Healing Abutment RP	Immediate Temporary Abutment RP	Procera Esthetic Abutment RP	Impression Coping Closed Tray RP
NobelSpeedy® Groovy RP					Plastic Coping Immediate Temporary Abutment RP	Esthetic Abutment RP	Impression Coping Open Tray RP
					QuickTemp™ Abutment Conical RP	15° Esthetic Abutment RP	Impression Coping Bar Brånemark System® RP
					Plastic Coping QuickTemp™ Abutment Conical RP	Snappy™ Abutment RP	Impression Coping Open Tray Multi-unit RP
					Temporary Abutment Engaging RP	Universal Base RP	Impression Coping Closed Tray Multi-unit RP
					GoldAdapt™ Engaging RP	Burn-out Coping RP	Impression Coping Closed Tray Multi-unit RP
					Temporary Abutment Non-Engaging RP	GoldAdapt™ Engaging RP	Impression Coping Bar Closed Tray Multi-unit RP
					Temporary Abutment Plastic Engaging RP	GoldAdapt™ Non-Engaging RP	Impression Coping Bar Closed Tray Multi-unit RP
					Temporary Abutment Plastic Non-Engaging RP	Multi-Unit Abutment RP	Temporary Coping Multi-unit Titanium RP
						17° Multi-Unit Abutment RP	Temporary Coping Multi-unit Plastic RP
						30° Multi-Unit Abutment RP	
						30° Multi-Unit Abutment Non-Engaging RP	
						Locator Abutment RP	
						Ball Abutment Titanium RP	
						Gold Abutment Bar Implant Level (vis clinique incluse) RP	
						Gold Coping Bar Multi-unit (vis de prothèse incluse) RP	
NobelSpeedy® Shorty WP	Cover Screw Brånemark System® WP	Implant Driver WP Brånemark System Implant Driver Wrench Adapter WP	Prosthetic Screw Multi-unit	Healing Abutment WP	Immediate Temporary Abutment WP	Procera Esthetic Abutment WP	Impression Coping Closed Tray WP
NobelSpeedy® Groovy WP					Plastic Coping Immediate Temporary Abutment WP	Esthetic Abutment WP	Impression Coping Open Tray WP
					QuickTemp™ Abutment Conical WP	15° Esthetic Abutment WP	Impression Coping Open Tray Multi-unit WP
					Plastic Coping QuickTemp™ Abutment Conical WP	Snappy™ Abutment WP	Impression Coping Closed Tray Multi-unit WP
					Temporary Abutment Engaging WP	Universal Base WP	Impression Coping Bar Closed Tray Multi-unit WP
					GoldAdapt™ Engaging WP	Burn-out Coping WP	Impression Coping Bar Closed Tray Multi-unit WP
					Temporary Abutment Non-Engaging WP	GoldAdapt™ Engaging WP	Temporary Coping Multi-unit Titanium WP
					Temporary Abutment Plastic Engaging WP	GoldAdapt™ Non-Engaging WP	Temporary Coping Multi-unit Plastic WP
					Temporary Abutment Plastic Non-Engaging WP	Multi-Unit Abutment WP	
						Locator Abutment WP	

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IFU1008. NobelReplace® Select™ TC

Nom du produit	Vis de couverture	Guide-implant	Vis cliniques	Piliers			Transfert d'empreinte	Autres
				de cicatrisation	provisoire	permanent		
Replace Select™ TC NP	Healing Screw RplSel TC NP Cover Screw NobRpl NP	Implant Driver NobRpl NP	Screw Multi-unit Angled Abutm NobRpl NP Screw Ceramic Abutment NobRpl NP Abutment Screw NobRpl NP	Healing Abutment NobRpl NP	Temporary Abutment Engaging NobRpl NP Temporary Abutment Non-Eng NobRpl NP Temporary Abutment Plastic Eng NobRpl NP Temp Abutm Plastic Non-Eng NobRpl NP Immediate Temporary Abutment NobRpl NP	Snappy™ Abutments NobRpl NP Esthetic Abutment NobRpl NP 15° Esthetic Abutment NobRpl NP Universal Base Tri-Channel NP GoldAdapt™ Engaging NobRpl NP GoldAdapt™ Non-Engaging NobRpl NP Multi-Unit Abutment NobRpl NP 17° Multi-unit Abutment NobRpl NP Ball Abutment Titanium NobRpl NP NobelProcera® Ab Ti NobelReplace® NP w-u NobelProcera® Ab Ti NobelReplace® NP Procera® Abutm Zr NobRpl NP	Impr Coping Open Tray NobRpl NP Ø 4,5 mm Impr Coping Closed Tray NobRpl NP Impr Cop Cl Tray Low Prof NRpl NP Impr Coping Open Tray NobRpl NP	PIB IBO Impl Retr Instr Hex & Tri-Ch NP/RP 22 mm Bone Mill with Guide NobRpl NP Ø 4,6 Implant Rescue Collar Tri-Channel Ø 3,5
Replace Select™ TC RP	Healing Screw RplSel TC RP Cover Screw NobRpl RP	Implant Driver NobRpl RP	Screw Multi-unit Angled Abutm NobRpl RP Screw Ceramic Abutment NobRpl RP/WP/6.0 Abutment Screw NobelReplace RP/WP/6.0	Healing Abutment NobRpl RP	Temporary Abutment Engaging NobRpl RP Temporary Abutment Non-Eng NobRpl RP Temporary Abutment Plastic Eng NobRpl RP Temp Abutm Plastic Non-Eng NobRpl RP Immediate Temporary Abutment NobRpl RP	Snappy™ Abutments NobRpl RP Esthetic Abutment NobRpl RP 15° Esthetic Abutment NobRpl RP Universal Base Tri-Channel RP GoldAdapt™ Engaging NobRpl RP GoldAdapt™ Non-Engaging NobRpl RP Multi-unit Abutment NobRpl RP 30° Multi-unit Abutment NobRpl RP 30° Multi-u Abutm Non-Eng NobRpl RP 17° Multi-unit Abutment NobRpl RP Ball Abutment Titanium NobRpl NP NobelProcera® Ab Ti NobelReplace® RP w-u NobelProcera® Ab Ti NobelReplace® RP Procera® Abutm Zr NobRpl RP	Impr Coping Closed Tray NobRpl RP Impr Cop Cl Tray Low Prof NRpl RP Impr Coping Open Tray NobRpl RP	PIB IBO Impl Retr Instr Hex & Tri-Ch NP/RP 31 mm Bone Mill with Guide NobRpl RP Ø 5,3 Implant Rescue Collar Tri-Channel Ø 4,3

IFU1010. NobelReplace® Conical Connection, NobelReplace® Conical Connection PMC, NobelReplace® Conical Connection TiUltra™

Nom du produit	Vis de couverture	Guide-implant	Vis cliniques	Piliers			Transfert d'empreinte	Autres
				de cicatrisation	provisoire	permanent		
NobelReplace® Conical Connection NP NobelReplace® Conical Connection PMC NP NobelReplace® Conical Connection TiUltra™ NP	Cover Screw CC NP	Implant driver CC NP	Vis clinique CC NP Omnigrip™ Clinical Screw CC NP On1™ Clinical Screw NP Screw MU Angled Abutment CC NP	Healing abutment CC NP Healing abutment bridge NP Slim healing abutment CC NP	Slim Temporary Abutment CC NP Temporary Snap Abutment CC NP Temporary Abutment Non Engaging CC NP	Universal Base CC NP Esthetic Abutment CC NP Procera Esthetic Abutment CC NP Multi-Unit Abutment Plus CC NP (0°/17°/30°) On1™ Base NP	Open Tray CC NP Closed Tray CC NP Bridge Open Tray CC NP	Guided Implant Mount NobRpl CC NP Bone Mill Guide CC NP Guide Pin CC NP NobelProcera® : NobelProcera® Abutment Zirconia CC NP NobelProcera® Abutment Ti CC NP NobelProcera® FCZ Implant Crown ASC CC NP Implant Bridge Restorations (Ti, Zr) Fixed et Fixed Removable Bars (Ti) Metal Adapter CC NP
NobelReplace® Conical Connection RP NobelReplace® Conical Connection PMC RP NobelReplace® Conical Connection TiUltra™ RP	Cover Screw CC RP	Implant driver CC RP	Vis clinique CC RP/WP Omnigrip™ Clinical Screw CC RP/WP On1™ Clinical Screw RP Screw MU Angled Abutment CC RP	Healing abutment CC RP Healing abutment bridge RP Slim healing abutment CC RP	Slim Temporary Abutment CC RP Temporary Snap Abutment CC RP Temporary Abutment Non Engaging CC RP	Universal Base CC NP Esthetic Abutment CC RP Procera Esthetic Abutment CC RP Multi-Unit Abutment Plus CC RP (0°/17°/30°) On1™ Base RP	Open Tray CC RP Closed Tray CC RP Bridge Open Tray CC RP	Guided Implant Mount NobRpl CC RP Bone Mill Guide CC RP Guide Pin CC RP NobelProcera® : NobelProcera® ASC Abutment Zirconia CC RP NobelProcera® Abutment Ti CC RP NobelProcera® FCZ Implant Crown ASC CC RP Implant Bridge Restorations (Ti, Zr) Fixed et Fixed Removable Bars (Ti) Metal Adapter CC RP

IFU1012. Replace Select™ Tapered TiUnite®, Replace Select™ Tapered PMC

Nom du produit	Vis de couverture	Guide-implant	Vis cliniques	Piliers			Transfert d'empreinte	Autres
				de cicatrisation	provisoire	permanent		
Replace Select™ Tapered NP	Cover Screw NobRpl NP Healing Screw RplSel TC NP	Implant Driver NobRpl NP	Abutment Screw NobelReplace® NP Screw Multi-unit Angled Abutment NobRpl NP Screw Ceramic Abutment NobRpl NP	Healing Abutment NobRpl NP	Temporary Abutment Engaging NobRpl NP Temporary Abutment Non-Eng NobRpl NP Immediate Temporary Abutment NobRpl NP Plastic/Temp Coping Esthetic Abutment NP Plastic/Temp Coping 15° Esthetic Ab NP	Snappy Abutment NobRpl NP Esthetic Abutment NobRpl NP 15° Esthetic Abutment NobRpl NP Universal Base Tri-Channel NP Universal Base Burn-out Coping NP/RP GoldAdapt™ Engaging NobRpl NP GoldAdapt™ Non-Engaging NobRpl NP Multi-Unit Abutment NobRpl NP 17° Multi-unit Abutment NobRpl NP Locator Abut for 3.5 RPL Sel Impl. Ball Abutment Titanium NobRpl NP NobelProcera® Ab Ti NobelReplace® NP Procera® Abutment Zr NobRpl NP	Impr Coping Open Tray NobRpl NP	PIB IBO Bone Mill with Guide NobRpl NP Bone Mill Guide NobelReplace® NP Impl Retr Instr Hex & Tri-Ch NP/RP Implant Rescue Collar Tri-Channel Ø 3,5
Replace Select™ Tapered RP	Cover Screw NobRpl RP Healing Screw RplSel TC RP	Implant Driver NobRpl RP	Abutment Screw NobelReplace RP/WP/6.0 Screw Multi-unit Angled Abutment NobRpl RP Screw Ceramic Abutment NobRpl RP/WP/6.0	Healing Abutment NobRpl RP	Temporary Abutment Engaging NobRpl RP Temporary Abutment Non-Eng NobRpl RP Immediate Temporary Abutment NobRpl RP	Snappy Abutment NobRpl RP Esthetic Abutment NobRpl RP 15° Esthetic Abutment NobRpl RP Universal Base Tri-Channel RP Universal Base Burn-out Coping NP/RP GoldAdapt™ Engaging NobRpl RP GoldAdapt™ Non-Engaging NobRpl RP Multi-unit Abutment NobRpl RP 17° Multi-unit Abutment NobRpl RP 30° Multi-unit Abutment NobRpl RP Locator Abut for 4.3 RPL Sel Impl. Ball Abutment Titanium NobRpl RP 3 mm NobelProcera® Ab Ti NobelReplace® RP Procera® Abutment Zr NobRpl RP	Impl Retr Instr Hex & Tri-Ch NP/RP	PIB IBO Bone Mill Guide NobelReplace® RP Bone Mill with Guide NobRpl RP Impl Retr Instr Hex & Tri-Ch NP/RP Implant Rescue Collar Tri-Channel Ø 4,3

Nom du produit	Vis de couverture	Guide-implant	Vis cliniques	Piliers			Transfert d'empreinte	Autres
				de cicatrisation	provisoire	permanent		
Replace Select™ Tapered WP	Cover Screw NobRpl WP	Implant Driver NobRpl WP	Screw Ceramic Abutment NobRpl RP/ WP/6.0	Healing Abutment NobRpl WP	Temporary Abutment Engaging NobRpl WP Temporary Abutment Non-Eng NobRpl WP Temp Abutment Plastic Non-Eng NobRpl WP Immediate Temporary Abutment NobRpl WP Plastic/Temp Cop Esthetic Ab NobRpl WP PI/Temp Coping 15° Esthetic Ab NobRpl WP	Snappy Abutment NobRpl WP Esthetic Abutment NobRpl WP 15° Esthetic Abutment NobRpl WP Universal Base Burn-out Coping WP Universal Base Tri-Channel WP GoldAdapt™ Non-Engaging NobRpl WP GoldAdapt™ Engaging NobRpl WP Multi-unit Abutment NobRpl WP Locator Abut for RPL Sel Impl. NobelProcera® Ab Ti NobelReplace® WP Procera® Abutment Zr NobRpl WP	Impr Coping Open Tray NobRpl WP	PIB IBO Bone Mill Guide NobelReplace® WP Bone Mill with Guide NobRpl WP Ø 6,5 Impl Retr Instr CC RP & Tri-Ch WP Impl Retr Instr CC RP& Tri-Ch WP& TCC RP
Replace Select™ Tapered 6.0	Cover Screw NobRpl 6,0 mm	Implant Driver NobRpl 6.0	Screw Ceramic Abutment NobRpl RP/ WP/6.0	Healing Abutment NobRpl 6.0	Temporary Abut Plastic Eng NobRpl 6.0 Temp Abutment Plastic Non-Eng NobRpl 6.0 Temporary Abutment Engaging NobRpl 6.0 Temporary Abutment Non-Eng NobRpl 6.0 Immediate Temporary Abutment NobRpl 6.0	Snappy Abutment NobRpl 6.0 Esthetic Abutment NobRpl 6.0 15° Esthetic Abutment NobRpl 6.0 GoldAdapt™ Non-Engaging NobRpl 6.0 GoldAdapt™ Engaging NobRpl 6.0 NobelProcera® Ab Ti NobelReplace® 6.0 Procera® Abutment Zr NobRpl 6.0	Impr Coping Open Tray NobRpl 6.0	PIB IBO Bone Mill Guide NobRpl 6.0 Bone Mill with Guide NobRpl 6.0 Ø 7 Impl Retr Instr Tri-Channel 6.0 22 mm

IFU1016. Vis de couverture

Vis de couverture	Type de connexion	Plate-forme implantaire	Système implantaire	Tournevis
Cover Screw Brånemark System® NP	Hexagonale externe	NP	Brånemark System® NobelSpeedy® Groovy	Cover Screw Driver Brånemark System Hexagon (IFU1085)
Cover Screw Brånemark System® RP		RP	Brånemark System® NobelSpeedy® Groovy NobelZygoma™ 0°	
Cover Screw Brånemark System® WP		WP	Brånemark System® NobelSpeedy® Groovy	
Cover Screw Conical Connection 3.0	Connexion conique	3.0	NobelActive®	Unigrip™ (IFU1085)
Cover Screw Conical Connection NP		NP	NobelActive® NobelParallel™ CC NobelReplace® CC	
Cover Screw Conical Connection RP		RP	NobelActive® NobelParallel™ CC NobelReplace® CC Trefoil™ systems	
Cover Screw Conical Connection WP		WP	NobelActive® NobelParallel™ CC NobelReplace® CC	
Cover Screw Nobel Biocare N1™ TCC NP	Connexion conique tri-ovale	NP	Nobel Biocare N1™ TiUltra™ TCC	Omnigrip™ Mini (IFU1085)
Cover Screw Nobel Biocare N1™ TCC RP		RP	Nobel Biocare N1™ TiUltra™ TCC	
Cover Screw NobelReplace® NP	Tri-rainure interne	NP	NobelReplace® Tapered Groovy, ReplaceSelect™	Unigrip™ (IFU1085)
Cover Screw NobelReplace® RP		RP		
Cover Screw NobelReplace® WP		WP		
Cover Screw NobelReplace® 6,0 mm		6,0		

IFU1018. Esthetic Abutment Conical Connection/Brånemark System®/ NobelReplace®, Esthetic Abutment Nobel Biocare N1™ TCC

Esthetic Abutment	Type de connexion	Vis clinique	Codage couleur	Tournevis	Chape en plastique/provisoire
Esthetic Abutment Conical Connection 3.0	Connexion conique	Clinical Screw Conical Connection 3.0	Aucun	Unigrip™	-
Esthetic Abutment Conical Connection NP	Connexion conique	Clinical Screw Conical Connection NP	●	Unigrip™	-
Esthetic Abutment Conical Connection RP	Connexion conique	Clinical Screw Conical Connection RP/ WP	●	Unigrip™	-
Esthetic Abutment Conical Connection WP	Connexion conique	Clinical Screw Conical Connection RP/ WP	●	Unigrip™	-
Esthetic Abutment Bmk System NP	Hexagonale externe	Abutment Screw Brånemark System® NP	Aucune	Unigrip™	Chape plastique/provisoire pour pilier Esthetic Abutment
Esthetic Abutment Bmk System RP	Hexagonale externe	Abutment Screw Brånemark System® RP	Aucune	Unigrip™	Chape plastique/provisoire pour pilier Esthetic Abutment
Esthetic Abutment Bmk System WP	Hexagonale externe	Abutment Screw Brånemark System® WP	Aucune	Unigrip™	Chape plastique/provisoire pour pilier Esthetic Abutment
Esthetic Abutment NobelReplace® NP	Tri-rainure	Abutment Screw NobelReplace® NP	●	Unigrip™	Chape plastique/provisoire pour pilier Esthetic Abutment
Esthetic Abutment NobelReplace® RP	Tri-rainure	Abutment Screw NobelReplace® RP/WP/6.0	●	Unigrip™	Chape plastique/provisoire pour pilier Esthetic Abutment

Esthetic Abutment	Type de connexion	Vis clinique	Codage couleur	Tournevis	Chape en plastique/provisoire
Esthetic Abutment NobelReplace® WP	Tri-rainure	Abutment Screw NobelReplace® RP/WP/6.0	●	Unigrip™	Chape plastique/provisoire pour pilier Esthetic Abutment
Esthetic Abutment NobelReplace® 6.0	Tri-rainure	Abutment Screw NobelReplace® RP/WP/6.0	●	Unigrip™	Chape plastique/provisoire pour pilier Esthetic Abutment
Esthetic Abutment Nobel Biocare N1™ TCC NP	Connexion conique tri-ovale	Clinical Screw Nobel Biocare N1 TCC NP	● (vis)	Omnigrip™ Mini	-
Esthetic Abutment Nobel Biocare N1™ TCC RP	Connexion conique tri-ovale	Clinical Screw Nobel Biocare N1 TCC RP	● (vis)	Omnigrip™ Mini	-

IFU1023. Piliers et bases Universal

Pilier/base Universal pour	Plates-formes disponibles	Hauteurs de limite cervicale disponibles	Couple de serrage	Tournevis
Connexion conique (CC)	NP	1,5 mm	35 Ncm	Unigrip™
	RP	3,0 mm		
	WP			
Connexion conique tri-ovale (TCC)	NP	1,5 mm	20 Ncm	Omnigrip™ mini
	RP	3,0 mm		
Tri-rainure	NP	1,5 mm	35 Ncm	Unigrip™
	RP	3,0 mm		
	WP			
Hexagonale externe	NP	1,5 mm	35 Ncm	Unigrip™
	RP	3,0 mm		
	WP			

IFU1040. Abutment Release Pins CC

Pilier	Abutment Release Pin CC 3.0	Abutment Release Pin CC NP	Abutment Release Pin CC RP/WP
Temporary Abutment Engaging Conical Connection	✓	s.o.	s.o.
Esthetic Abutment Conical Connection	✓	s.o.	s/o
GoldAdapt™ Engaging Conical Connection	✓	s.o.	s/o
NobelProcera® Abutment Zirconia Conical Connection	s/o	✓	✓

IFU1043. Abutment Screw Retrieval Instrumentation

Instrument	Type de connexion	Plate-forme	Handle for Machine instruments/ Prosthetic Torque Wrench Adapter	Contre-angle (compatible avec DIN EN ISO 1797)
Abutment Screw Remove 3.0	Conical Connection External Hex Internal Tri-Channel	3.0	-	X
		NP	-	X
Abutment Screw Remove NP	Trioval Conical Connection (NP and RP platform only)	RP	-	X
		WP	-	X
		6.0	-	X
		3.0	-	X
Abutment Screw Remove RP/WP/6.0	Conical Connection External Hex Internal Tri-Channel	NP	-	X
		RP	-	X
		WP	-	X
Abutment Screw Retrieval Reverse Drill 3.0/NP	Conical Connection External Hex Internal Tri-Channel	6.0	-	X
		RP	-	X
		WP	-	X
Abutment Screw Retrieval Reverse Drill RP/WP/6.0	Conical Connection External Hex Internal Tri-Channel	6.0	-	X
		NP	-	X
		RP	-	X
		WP	-	X
Abutment Screw Retrieval Reverse Drill CC 3.0/NP & TCC NP/RP	Conical Connection External Hex Internal Tri-Channel	6.0	-	X
		NP	-	X
		RP	-	X
		WP	-	X
Rescue Drill Guide Conical Connection 3.0	Connexion conique tri-ovale	NP	-	X
		RP	-	X
Rescue Drill Guide Conical Connection NP	Connexion conique	3.0	-	-
		NP	-	-
Rescue Drill Guide Conical Connection RP	Connexion conique	RP	-	-
		WP	-	-
Rescue Drill Guide Conical Connection WP	Hexagonale externe	NP	-	-
		RP	-	-
Rescue Drill Guide Hexagonale externe NP	Hexagonale externe	RP	-	-
		WP	-	-
Rescue Drill Guide Hexagonale externe RP	Hexagonale externe	NP	-	-
		WP	-	-
Rescue Drill Guide Hexagonale externe WP	Hexagonale externe	RP	-	-
		WP	-	-
Rescue Drill Guide Tri-rainure interne Tri-Channel NP	Tri-rainure interne	NP	-	-
		RP	-	-
		WP	-	-
Rescue Drill Guide Tri-rainure RP	Tri-rainure interne	NP	-	-
		WP	-	-
Rescue Drill Guide Tri-rainure WP	Tri-rainure interne	6.0	-	-
		NP	-	-
Rescue Drill Guide Tri-rainure 6.0	Tri-rainure interne	RP	-	-
		WP	-	-
Rescue Drill Guide Nobel Biocare N1™ TCC NP	Connexion conique tri-ovale	NP	-	-
		RP	-	-
Rescue Drill Guide Nobel Biocare N1™ TCC RP	Connexion conique tri-ovale	RP	-	-
		WP	-	-

Instrument	Type de connexion	Plate-forme	Handle for Machine instruments/ Prosthetic Torque Wrench Adapter	Contre-angle (compatible avec DIN EN ISO 1797)
Abutment Screw Retrieval Instrument 3.0/NP	Conical Connection External Hex Internal Tri-Channel	3.0	X	X
		NP	X	X
		RP	X	X
Abutment Screw Retrieval Instrument RP/WP/6.0		WP	X	X
		6.0	X	X
Abutment Screw Retrieval Instrument CC 3.0/NP & TCC NP/RP	Conical Connection External Hex Internal Tri-Channel	3.0	X	X
		NP	X	X
		RP	X	X
		WP	X	X
	6.0	X	X	
	Connexion conique tri-ovale	NP	X	X
		RP	X	X
Screw Tap Repair M1.4	Connexion conique	3.0	-	X
Screw Tap Repair M1.6	Connexion conique Hexagonale externe	NP	-	X
Screw Tap Repair M1.8	Tri-rainure interne	NP	-	X
Screw Tap Repair M2.0	Connexion conique	RP	-	X
		WP	-	X
	Tri-rainure interne	RP	-	X
		WP	-	X
		6.0	-	X
	Hexagonale externe	RP	-	X
Screw Tap Repair M2.5	Hexagonale externe	WP	-	X
Screw Tap Repair Tool Nobel Biocare N1™ TCC NP	Connexion conique tri-ovale	NP	-	X
Screw Tap Repair Tool Nobel Biocare N1™ TCC RP		RP	-	X
Nobel Biocare N1™ Base Screw Removal Tool NP/RP		NP/RP	X	-

IFU1048. Forets anatomiques Tapered Drill, forets pour os de forte densité Dense Bone Drills, tarauds Screw Taps

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IFU1074. Concept On1™

Nom du produit	Implant	Pilier	Vis	Tournevis
On1™ Base CC NP 1,75 mm	Implants avec connexion conique interne : - NobelActive® - NobelActive® TiUltra™ - NobelReplace® CC - NobelReplace® CC TiUltra™ - NobelParallel™ CC - NobelParallel™ CC TiUltra™	On1™ Temporary Abutment Engaging On1™ Temporary Abutment Non-Engaging On1™ Universal Abutment On1™ Universal Abutment Non-Engaging On1™ Esthetic Abutment	On1™ Clinical Screw 1,75 mm NP	On1™ Screwdriver
On1™ Base CC NP 2,5 mm			On1™ Clinical Screw 2,5 mm NP	
On1™ Base CC RP 1,75 mm			On1™ Clinical Screw 1,75 mm RP	
On1™ Base CC RP 2,5 mm			On1™ Clinical Screw 2,5 mm RP	
On1™ Base CC WP 1,75 mm			On1™ Clinical Screw 1,75 mm WP	
On1™ Base CC WP 2,5 mm			On1™ Clinical Screw 2,5 mm WP	
On1™ Universal Abutment 0,3 mm NP	s.o.	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw	Unigrip™ Screwdriver
On1™ Universal Abutment 1,25 mm NP			On1™ Prosthetic Lab Screw	
On1™ Universal Abutment 0,3 mm RP				
On1™ Universal Abutment 1,25 mm RP				
On1™ Universal Abutment 0,3 mm WP				
On1™ Universal Abutment 1,25 mm WP				
On1™ Esthetic Abutment Titanium 0,3 mm NP	s.o.	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw	Unigrip™ Screwdriver
On1™ Esthetic Abutment Titanium 0,3 mm RP			On1™ Prosthetic Lab Screw	
On1™ Esthetic Abutment Titanium 0,3 mm WP				
On1™ Base Xeal™ CC NP 1,75 mm	Implants avec connexion conique interne : - NobelActive® - NobelActive® TiUltra™ - NobelReplace® CC - NobelReplace® CC TiUltra™ - NobelParallel™ CC - NobelParallel™ CC TiUltra™	On1™ Temporary Abutment Engaging On1™ Temporary Abutment Non-Engaging On1™ Universal Abutment On1™ Universal Abutment Non-Engaging On1™ Esthetic Abutment	On1™ Clinical Screw 1,75 mm NP	On1™ Screwdriver
On1™ Base Xeal™ CC NP 2,5 mm			On1™ Clinical Screw 2,5 mm NP	
On1™ Base Xeal™ CC RP 1,75 mm			On1™ Clinical Screw 1,75 mm RP	
On1™ Base Xeal™ CC RP 2,5 mm			On1™ Clinical Screw 2,5 mm RP	
On1™ Base Xeal™ CC WP 1,75 mm			On1™ Clinical Screw 1,75 mm WP	
On1™ Base Xeal™ CC WP 2,5 mm			On1™ Clinical Screw 2,5 mm WP	

Nom du produit	Implant	Pilier	Vis	Tournevis
On1™ Universal Abutment Non-Eng NP 0,3 mm	s.o.	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw On1™ Prosthetic Lab Screw	Unigrip™ Screwdriver
On1™ Universal Abutment Non-Eng NP 1,25 mm				
On1™ Universal Abutment Non-Eng RP 0,3 mm				
On1™ Universal Abutment Non-Eng RP 1,25 mm				
On1™ Universal Abutment Non-Eng WP 0,3 mm				
On1™ Universal Abutment Non-Eng WP 1,25 mm				
On1™ Healing Cap NP 1,5 mm	s.o.	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw	Unigrip™ Screwdriver
On1™ Healing Cap NP 2,5 mm				
On1™ Healing Cap RP 1,5 mm				
On1™ Healing Cap RP 2,5 mm				
On1™ Healing Cap WP 1,5 mm				
On1™ Healing Cap WP 2,5 mm				
On1™ IOS Healing Cap 6 mm NP	s.o.	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw	Unigrip™ Screwdriver
On1™ IOS Healing Cap 4,5 mm NP				
On1™ IOS Healing Cap 6 mm RP				
On1™ IOS Healing Cap 4,5 mm RP				
On1™ IOS Healing Cap 5 mm WP				
On1™ IOS Healing Cap 4 mm WP				
On1™ Temporary Abutment Engaging NP	s.o.	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw On1™ Prosthetic Lab Screw	Unigrip™ Screwdriver
On1™ Temporary Abutment Non-Engaging NP				
On1™ Temporary Abutment Engaging RP				
On1™ Temporary Abutment Non-Engaging RP				
On1™ Temporary Abutment Engaging WP				
On1™ Temporary Abutment Non-Engaging WP				
On1™ Impression Coping Closed Tray NP	s.o.	On1™ Base On1™ Base Xeal	s.o.	Unigrip™ Screwdriver
On1™ Impression Coping Closed Tray RP				
On1™ Impression Coping Closed Tray WP				
On1™ Impression Coping Open Tray NP				
On1™ Impression Coping Open Tray RP				
On1™ Impression Coping Open Tray WP				
On1™ Impression Coping Open Tray N-Eng NP	s.o.	On1™ Base On1™ Base Xeal	s.o.	Unigrip™ Screwdriver
On1™ Impression Coping Open Tray N-Eng RP				
On1™ Impression Coping Open Tray N-Eng WP				
On1™ Clinical Screw 1,75 mm NP	Implants avec connexion conique interne : - NobelActive® - NobelActive® TiUltra™ - NobelReplace® CC - NobelReplace® CC TiUltra™ - NobelParallel™ CC - NobelParallel™ CC TiUltra™	On1™ Temporary Abutment Engaging On1™ Temporary Abutment Non-Engaging On1™ Universal Abutment On1™ Universal Abutment Non-Engaging On1™ Esthetic Abutment	On1™ Prosthetic Screw On1™ Prosthetic Lab Screw On1™ Impression Coping Screw	On1™ Screwdriver
On1™ Clinical Screw 2,5 mm NP				
On1™ Clinical Screw 1,75 mm RP				
On1™ Clinical Screw 2,5 mm RP				
On1™ Clinical Screw 1,75 mm WP				
On1™ Clinical Screw 2,5 mm WP				
On1™ Prosthetic Screw	s.o.		On1™ Clinical Screw	Unigrip™ Screwdriver
On1™ Screwdriver Machine 21 mm	s.o.	On1™ Base On1™ Base Xeal	On1™ Clinical Screw	s.o.
On1™ Screwdriver Manual 25 mm	s.o.	On1™ Base On1™ Base Xeal	On1™ Clinical Screw	s.o.

IFU1075. Piliers Multi-unit Abutments et capuchons de cicatrisation Healing Caps

Piliers Multi-unit Abutment/ plates-formes	Système implantaire	Capuchon de cicatrisation	Tournevis
Connexion conique interne (CC)			
Multi-unit Abutment Plus CC NP/RP/WP Multi-unit Abutment Xeal™ CC NP/RP/WP 17° Multi-unit Abutment Plus CC NP/RP/WP 17° Multi-unit Abutment Xeal™ CC NP/RP/WP 30° Multi-unit Abutment Plus CC NP/RP 30° Multi-unit Abutment Xeal™ CC NP/RP	NobelActive® NobelParallel™ CC NobelReplace™ CC	Healing Cap Multi-unit Titanium	Capuchon de cicatrisation : Unigrip™ MUA droits : MUA Driver MUA angulés : Unigrip™
Connexion conique tri-ovale (TCC)			
Multi-unit Abutment Xeal™ Nobel Biocare N1™ TCC NP/RP 17° Multi-unit Abutment Xeal™ Nobel Biocare N1™ TCC NP/RP 30° Multi-unit Abutment Xeal™ Nobel Biocare N1™ TCC RP	Nobel Biocare N1™	Healing Cap Multi-unit Titanium	Capuchon de cicatrisation : Unigrip™ MUA droits : MUA Driver MUA angulés : Omnigrip™ Mini
Connexion tri-rainure interne			
Multi-Unit Abutment NobRpl NP/RP/WP 17° Multi-Unit Abutment NobelReplace® NP/RP 30° Multi-Unit Abutment NobelReplace® RP	Replace Select™ NobelSpeedy® Replace™ NobelSpeedy® Replace	Healing Cap Multi-unit Titanium	Capuchon de cicatrisation : Unigrip™ MUA angulés : Unigrip™
Hexagonale externe			
Multi-unit Abutments Brånemark System® NP/RP 17° Multi-Unit Abutments Brånemark System® NP/RP 30° Multi-Unit Abutments Brånemark System® RP	Brånemark System® NobelSpeedy® Groovy	Healing Cap Multi-unit Titanium	Capuchon de cicatrisation : Unigrip MUA droits : MUA Driver MUA angulés : Unigrip
Multi-unit Abutments Brånemark System® WP	Brånemark System® NobelSpeedy® Groovy	Capuchon de cicatrisation Multi-unit Brånemark System® WP	Capuchon de cicatrisation : Unigrip MUA droits : MUA Driver
Multi-unit Abutments Brånemark System® Zygoma 17° Multi-unit Abutments Brånemark System® Zygoma	Brånemark System® Zygoma	Healing Cap Multi-unit Titanium	Capuchon de cicatrisation : Unigrip MUA droits : MUA Driver MUA angulés : Unigrip
45° Multi-unit Abutments External Hex RP 60° Multi-unit Abutments External Hex RP	Implant NobelZygoma™ 0°		Capuchon de cicatrisation : Unigrip MUA angulés : Unigrip

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Handle for Implant Rescue Collar and Drill Guides

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Implant Driver Wrench Adapters Brånemark System® NP/RP/WP

Implant	Instructions d'utilisation	Clé à torque
Brånemark System®	IFU1007	Brånemark System® Manual Torque Wrench
NobelSpeedy®	IFU1007	

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IFU1096. Abutment Retrieval Instrumentation

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Clé	Adaptateur pour clé	Instruments
NobelActive® Manual Torque Wrench Surgical	NOBELREPLACE® Manual Torque Wrench Adapter Surgical	Implant Driver CC 3.0 28 mm (IFU1090)
		Implant Driver CC 3.0 pour Slim Abutment Implant Driver CC NP (IFU1090) Implant Driver CC NP pour Slim Abutment (IFU1090)
NobelReplace® Manual Torque Wrench Surgical	NOBELREPLACE® Manual Torque Wrench Adapter Surgical	Implant Driver CC RP (IFU1090)
		Implant Driver CC RP pour Slim Abutment (IFU1090) Implant Driver CC WP (IFU1090) Implant Retrieval Instruments (IFU1097) Abutment Screw Retrieval Instrument (IFU1043)
Brånemark System® Manual Torque Wrench Surgical	Brånemark System® Manual Torque Wrench Surgical	Implant Driver Brånemark Syst NP (IFU1090)
		Implant Driver Brånemark Syst RP (IFU1090) Implant Driver Brånemark Syst WP (IFU1090)
Manual Torque Wrench Surgical Nobel Biocare N1™	NobelReplace® Manual Torque Wrench Adapter Surgical	Implant Driver Nobel Biocare N1™ TCC NP (IFU1087)
		Implant Driver Nobel Biocare N1™ TCC RP (IFU1087) Implant Retrieval Instruments (IFU1097) Abutment Screw Retrieval Instrument (IFU1043)
Manual Torque Wrench Prosthetic	Manual Torque Wrench Adapter Prosthetic	Screwdriver Machine Unigrip (IFU1085)
		Screwdriver Machine Multi-Unit (IFU1085) Omnigrip™ Screwdriver Machine (IFU1085) On1™ Base Screwdriver Machine (IFU1074) Screwdriver Machine Multi-unit Brånemark System (IFU1085) Screwdriver Machine Ball Abutment (IFU1085)
Manual Torque Wrench Prosthetic Nobel Biocare N1™	Manual Torque Wrench Adapter Prosthetic	Screwdriver Machine Omnigrip™ Mini (IFU1085)
		Screwdriver Machine Multi-Unit (IFU1085)

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Pilier original	Tournevis de rechange
Pilier boule	Screwdriver Machine Ball Abutment Screwdriver Manual Ball Abutment

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Nom du produit	Implant	Sleeve	Driver	Autres
Guided Surgery Implant Mount NobelActive® NP	NobelActive®/ NobelActive® TiUltra™ NP	Guided Sleeve NP	Surgical Driver	Manual Torque Wrench Adapter Surgical
Guided Surgery Implant Mount NobelActive® RP 4.3	NobelActive®/ NobelActive® TiUltra™ RP 4.3	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	Connection to Handpiece
Guided Surgery Implant Mount NobelActive® RP 5.0	NobelActive®/ NobelActive® TiUltra™ RP 5.0	Guided Sleeve 6.0/WP		
Guided Surgery Implant Mount NobelActive® WP 5.5	NobelActive®/ NobelActive® TiUltra™ WP 5.5			
Guided Implant Mount Bmk Syst NP	Brånemark System® NP	Guided Sleeve NP	Surgical Driver	Manual Torque Wrench Adapter Surgical Brånemark System®
Guided Implant Mount Bmk Syst RP	Brånemark System® RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	Connection to handpiece
Guided Implant Mount Bmk Syst WP	Brånemark System® WP	Guided Sleeve 6.0/WP		
Guided Implant Mount NobRpl NP	NobelReplace® NP	Guided Sleeve NP	Surgical Driver	NobelReplace® Manual Torque Wrench Adapter Surgical Connection to handpiece
Guided Implant Mount NobRpl RP	NobelReplace® RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	
Guided Implant Mount NobRpl WP	NobelReplace® WP	Guided Sleeve 6.0/WP	Guided Sleeve 6.0/WP	
Guided Implant Mount NobRpl 6.0	NobelReplace® 6.0	Guided Sleeve 6.0/WP	Guided Sleeve 6.0/WP	
Guided Implant Mount NobelReplace® Conical Connection NP 3.5	NobelReplace® CC NP	Guided Sleeve NP	Surgical Driver	NobelReplace® Manual Torque Wrench Adapter Surgical Connection to handpiece
Guided Implant Mount NobelReplace® Conical Connection RP 4.3	NobelReplace® CC RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	
Guided Implant Mount NobelReplace® Conical Connection RP 5.0	NobelReplace® CC RP 5.0	Guided Sleeve 6.0/WP		
Guided Implant Mount NobelParallel™ CC NP	NobelParallel™ CC NP	Guided Sleeve NP	Surgical Driver	NobelReplace® Manual Torque Wrench Adapter Surgical Connection to handpiece
Guided Implant Mount NobelParallel™ CC RP 4.3	NobelParallel™ CC RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	
Guided Implant Mount NobelParallel™ CC RP 5.0	NobelParallel™ CC RP 5.0	Guided Sleeve 6.0/WP		
Guided Implant Mount NobelParallel™ CC WP 5.5	NobelParallel™ CC WP	Guided Sleeve 6.0/WP		

Nom du produit	Implant	Sleeve	Driver	Autres
Guided Template Abutm w Screw NobelActive® NP	NobelActive®/ NobelActive® TiUltra™ NP	Guided Sleeve NP	Screwdriver Manual/ Machine Unigrip™	s.o.
Guided Template Abutm w Screw NobelActive® RP 4.3	NobelActive®/ NobelActive® TiUltra™ RP 4.3	Guided Sleeve RP		
Guided Template Abutm w Screw NobelActive® RP 5.0	NobelActive®/ NobelActive® TiUltra™ RP 5.0	Guided Sleeve 6.0/WP		
Guided Template Abutment with Screw NobRpl 6.0	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ 6.0 et Replace Select™ Tapered TiUnite® 6.0			
Guided Template Abutment with Screw NobRpl NP	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ NP et Replace Select™ Tapered TiUnite® NP	Guided Sleeve NP		
Guided Template Abutment with Screw NobRpl RP	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ RP et Replace Select™ Tapered TiUnite® RP	Guided Sleeve RP		
Guided Template Abutment with Screw NobRpl WP	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ WP et Replace Select™ Tapered TiUnite® WP	Guided Sleeve 6.0/WP		
Guided Template Abutment CC NP	Implants NobelActive®/ NobelActive® TiUltra™ NP, NobelParallel™ CC/ NobelParallel™ CC TiUltra™ NP, NobelReplace® CC/ NobelReplace® CC TiUltra™ NP, NobelReplace® CC PMC NP	Guided Sleeve NP		
Guided Template Abutment CC NP	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ et CC PMC NP			
Guided Template Abutment CC RP 4.3	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ & CC PMC RP 4.3	Guided Sleeve RP		
Guided Template Abutment CC RP 5.0	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ & CC PMC RP 5.0	Guided Sleeve 6.0/WP		
Guided Template Abutment CC WP 5.5	Implants NobelReplace® CC/ NobelReplace® CC TiUltra™ & CC PMC WP 5.5			
Guided Template Abutment with Screw Brånemark System® NP	Implants Brånemark System® MK III NP et NobelSpeedy® Groovy NP	Guided Sleeve NP		
Guided Template Abutment avec Screw Brånemark System® RP	Implants Brånemark System® MK III RP et NobelSpeedy® Groovy RP	Guided Sleeve RP		
Guided Template Abutment with Screw Brånemark System® WP	Implants Brånemark System® MK III WP et NobelSpeedy® Groovy WP	Guided Sleeve 6.0/WP		
Guided Template Abutment with Screw NobelActive® WP 5.5	NobelActive®/ NobelActive® TiUltra™ WP 5.5			
Handle for Guided Drill Guide	s.o.	s.o.	s.o.	Guided Drill Guides

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Nom du produit	Connexions des contre-angles
Guided Start Drill/Counterb Bmk Syst NP	EN ISO 1797
Guided Start Drill/Counterbore NobRpl NP	
Guided Start Drill/Counterb Bmk Syst RP	
Guided Start Drill/Counterbore NobRpl RP	
Guided Start Drill/Counterb BmkSyst WP Ø 5	
Guided Start Drill/Counterbore NobRpl WP	
Guided Start Drill/Counterb BmkSyst WP Ø 6	
Guided Start Drill/Counterb NobRpl 6.0	
Guided Twist Drill 2x(10+)7-18 mm	
Guided Twist Drill 2,8x(10+)7-18 mm	
Guided Twist Drill 3x(10+)7-18 mm	
Guided Twist Drill 3,2x(10+)7-18 mm	
Guided Twist Drill 3,4x(10+)7-18 mm	
Guided Twist Drill 3,8x(10+)7-18 mm	
Guided Twist Drill 4,2x(10+)7-18 mm	
Guided Twist Drill 5x(10+)7-18 mm	
Guided Dense Bone Drill Taper NP 3,5x13 mm	
Guided Dense Bone Drill Taper NP 3,5x16 mm	
Guided Dense Bone Drill Taper RP 4,3x13 mm	
Guided Dense Bone Drill Taper RP 4,3x16 mm	
Guided Dense Bone Drill Taper WP 5x13 mm	
Guided Dense Bone Drill Taper WP 5x16 mm	
Guided Dense Bone Drill Taper 6,0 6x13 mm	
Guided Dense Bone Drill Taper 6,0 6x16 mm	
Guided Screw Tap NP 10-13 mm	
Guided Screw Tap RP Ø 3,75 7-13 mm	
Guided Screw Tap Tapered RP	
Guided Screw Tap WP Ø 5 7-13 mm	
Guided Screw Tap Tapered WP	
Guided Screw Tap 6.0/WP Ø 6 7-13 mm	
Guided Screw Tap Tapered 6.0	
Guided StartDrill/Countb former MkIII RP	
Guided Screw Tap RP Ø 4 7-13 mm	
Guided Twist Drill 3,2x(10+)7-13 mm	
Guided Twist Drill 3,4x(10+)7-13 mm	
Guided Twist Drill 3,8x(10+)7-13 mm	
Guided Twist Drill 4,2x(10+)7-13 mm	
Guided Twist Drill 5x(10+)7-13 mm	
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 2,4/2,8 (10+) 7-13 mm	
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 2,8/3,2 (10+) 7-13 mm	
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 3,2/3,6 (10+) 7-13 mm	
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 3,8/4,2 (10+) 7-13 mm	
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 4,2/4,6 (10+) 7-13 mm	
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 2,4/2,8 (10+) 7-18 mm	
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 2,8/3,2 (10+) 7-18 mm	

Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 3,2/3,6 (10+) 7-18 mm
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 3,8/4,2 (10+) 7-18 mm
Guided Screw Tap NobelActive® NP 8,5-10 mm
Guided Screw Tap NobelActive® RP 4.3 8,5-10 mm
Guided Screw Tap NobelActive® RP 5.0 8,5-10 mm
Guided Twist Drill Tapered 2 x (10+) 8-16 mm
Foret de démarrage guidé Guided Start Drill
Guided Counterbore NobelReplace® 3.5
Guided Counterbore NobelReplace® 4.3
Guided Counterbore NobelReplace® 5.0
Guided Counterbore NobelReplace® 6.0
Guided Screw Tap Tapered NP
Guided Pilot Twist Drill 1,5x(10+)10-15 mm
Foret hélicoïdal à paliers guidé Guided Twist Step Drill 4,2/5,0 (10+) 7-13 mm
Foret hélicoïdal à paliers guidé Guided Twist Step Drill 4,2/5,0 (10+) 7-18 mm
Guided Screw Tap NobelActive® WP 5.5 7-8.5 mm
Guided Screw Tap NobelActive® WP 5.5 10-15 mm
Foret hélicoïdal à paliers guidé Guided Twist Step Drill Ø 4,2/4,6x(10+)7-18 mm
Taraud guidé NobelParallel™ CC 3.75 7-10 mm
Taraud guidé NobelParallel™ CC 3.75 11,5-13 mm
Taraud guidé NobelParallel™ CC 3.75 15-18 mm
Taraud guidé NobelParallel™ CC 4.3 7-10 mm
Taraud guidé NobelParallel™ CC 4.3 11,5-13 mm
Taraud guidé NobelParallel™ CC 4.3 15-18 mm
Taraud guidé NobelParallel™ CC 5.0 7-10 mm
Taraud guidé NobelParallel™ CC 5.0 11,5-13 mm
Taraud guidé NobelParallel™ CC 5.0 15-18 mm
Taraud guidé NobelParallel™ CC 5.5 7-10 mm
Taraud guidé NobelParallel™ CC 5.5 11,5-13 mm
Taraud guidé NobelParallel™ CC 5.5 15-18 mm
Guided Tissue Punch NP-sleeve (3.75)
Guided Tissue Punch WP-sleeve (5.5)
Guided Counterbore NobelParallel™ CC 3.75
Guided Counterbore NobelParallel™ CC 4.3
Guided Counterbore NobelParallel™ CC 5.0
Guided Counterbore NobelParallel™ CC 5.5
Guided Dense Bone Screw Tap NobelActive® NP 8,5-18 mm
Guided Dense Bone Screw Tap NobelActive® RP 4.3 8,5-18 mm
Guided Dense Bone Screw Tap NobelActive® RP 5.0 8,5-18 mm
Guided Dense Bone Screw Tap NobelActive® WP 5.5 7-8,5 mm
Guided Dense Bone Screw Tap NobelActive® WP 5.5 10-15 mm
Guided Drill Tapered NP
Guided Drill Tapered RP
Guided Drill Tapered WP
Guided Drill Tapered 6.0

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IFU3018. Bridge implantaire NobelProcera® en zircone

	Type de connexion											
	Hexagonale externe			Tri-rainure interne				CC			MUA	
	NP	RP	WP	NP	RP	WP	6.0	NP	RP	WP	NP/RP	WP
Bridge implantaire NobelProcera en zircone												
Adaptateur clinique	s.o.							Adaptateur en métal pour implant Bridge Clinical CC			s.o.	
Adaptateur de laboratoire	s.o.							Adaptateur en métal pour implant Bridge Lab CC			s.o.	
Vis clinique	Screw Ceramic Abutment Brånemark System®							Omnigrip™ Clinical Screw CC			Prosthetic Screw Multi-unit Abutment Omnigrip™ Mini	
Vis de laboratoire	Lab Screw Implant Level External Hex			s.o. (*)				Omnigrip™ Lab Screw CC Lab Screw Implant Level Conical Connection			Lab Screw Multi-unit Abutment Omnigrip™ Mini	
Implant/répliques MUA	Implant Replica Brånemark System® Elos Accurate Analog pour modèles imprimés Hexagonale externe			Implant Replicas NobelReplace® Elos Accurate Analog pour modèles imprimés Tri-rainure				Implant Replica Conical Connection IOS Model Implant Replica CC			Multi-unit Abutment Plus Replica Elos Accurate Analog pour modèles imprimés Multi-unit	Multi-unit Abutment Plus Replica Abutment Replica Multi-unit Bmk Syst Elos Accurate Analog pour modèles imprimés Multi-unit
Analogues de protection	Protection Analog Bmk Syst			Protection Analog NobelReplace®				Protection Analog/Drill Guide CC			Protection Analog Multi-unit	Protection Analog Multi-unit Protection Analog M-u Bmk Syst
MUA	Multi-unit Abutment Bmk Syst 17° Multi-unit Abutment Bmk Syst 30° Multi-unit Abutment Bmk Syst Zygoma Abutment Brånemark System® Zygoma Zygoma 17° Brånemark System® Zygoma 17° 45° Multi-unit Abutment Ext Hex 60° Multi-unit Abutment Ext Hex			Multi-unit Abutment NobRpl 17° Multi-unit Abutment NobRpl 30° Multi-unit Abutment NobRpl				Multi-unit Abutment Plus Conical Connection 17° Multi-unit Abutment Plus Conical Connection 30° Multi-unit Abutment Plus Conical Connection Multi-unit Abutment Xeal™ Conical Connection 17° Multi-unit Abutment Xeal™ Conical Connection 30° Multi-unit Abutment Xeal™ Conical Connection			s.o.	
Implants	NobelSpeedy® Groovy NobelSpeedy® Shorty	Brånemark System® Mk III Brånemark System® Mk III TiUnite® NobelSpeedy® Groovy NobelSpeedy® Shorty Zygoma Implant NobelZygoma™	NobelSpeedy® Groovy NobelSpeedy® Shorty	Replace Select™ TC				NobelActive® NobelActive® TiUltra™ NobelParallel™ CC NobelParallel™ CC TiUltra™ NobelReplace® CC NobelReplace® CC PMC NobelReplace® CC TiUltra™	NobelActive® NobelActive® TiUltra™ NobelParallel™ CC NobelParallel™ CC TiUltra™ NobelReplace® CC NobelReplace® CC PMC NobelReplace® CC TiUltra™	NobelActive® NobelActive® TiUltra™ NobelParallel™ CC NobelParallel™ CC TiUltra™	s.o.	

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