

# Healing Abutments



## Important – Disclaimer of Liability

This product is part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Nobel Biocare. Non-recommended use of products made by third parties in conjunction with Nobel Biocare products will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare products has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The user is also obliged to study the latest developments in regard to this Nobel Biocare product and its applications regularly. In cases of doubt, the user has to contact Nobel Biocare. Since the utilization of this product is under the control of the user, they are his/her responsibility. Nobel Biocare does not assume any liability whatsoever for damage arising thereof.

Please note that some products detailed in this Instruction for Use may not be regulatory cleared, released or licensed for sale in all markets.

## Description

Healing abutments are premanufactured dental implant abutments which can be directly connected to the endosseous dental implant to support healing of the surrounding soft tissue.

## Healing Abutments / Healing Screw

- Healing Abutments Conical Connection are available in NP/RP/WP platforms, feature a conical connection and can be used with Nobel Biocare's NobelActive®, NobelParallel™ CC and/or NobelReplace® CC implant systems.
- Healing Abutments Bridge Conical Connection are available in NP/RP/WP platforms, feature a conical connection and can be used with Nobel Biocare's NobelActive®, NobelParallel™ CC and/or NobelReplace® CC implant systems.
- Healing Abutments Nobel Biocare N1™ TCC are available in NP/RP platforms and feature a tri-oval conical connection and can be used with Nobel Biocare's Nobel Biocare N1™ implant system.
- Healing Abutments NobelReplace® are available in NP/RP/WP and 6.0 platform, feature an internal tri-channel connection and can be used with Nobel Biocare's NobelReplace®, Replace Select™ and/or NobelSpeedy® Replace™ implant systems.
- Healing Screw Replace Select™ is available in NP and RP platforms, feature an internal tri-channel connection and can be used with Replace Select™ TC.
- Healing Abutments Brånemark System® are available in NP/RP/WP platforms, feature an external hex connection and can be used with Nobel Biocare's Brånemark System® and/or NobelSpeedy® Groovy® implant systems. The RP can be used with NobelZygoma.
- Healing Abutments Brånemark System® Zygoma feature an external hex connection and can be used with Nobel Biocare's Brånemark System® Zygoma implant system.

## Slim Healing Abutments

- Slim Healing Abutments Conical Connection are available in 3.0/NP/RP platforms, feature a conical connection and can be used with Nobel Biocare's NobelActive®, NobelParallel™ CC and/or NobelReplace CC implant systems.

## Healing Abutments Anatomical PEEK

- Healing Abutments Anatomical PEEK Conical Connection are available in WP platform, feature a conical connection and can be used with Nobel Biocare's NobelActive® and/or NobelParallel™ CC implant systems.

The following tables summarize the implant platforms which are compatible with the various healing abutments, including the specifications for required screwdrivers, and other key information for each type of healing abutment, based on their connection type.

Healing Abutment for	Available platforms	Color coding	Screwdriver
Conical Connection (CC)	3.0	None	Unigrip™
	NP	●	
	RP	●	
	WP	●	
Tri-oval Conical Connection (TCC)	NP	● (screw)	OmniGrip™ Mini
	RP	● (screw)	
Tri-channel	NP	●	Unigrip™
	RP	●	
	WP	●	
	6.0	●	
External Hex	NP	None	Unigrip™
	RP	None	
	WP	None	

**Table 1 – Healing Abutments – Compatible Implant Platforms and Screwdrivers**

Healing Screw for	Available platforms	Color coding	Screwdriver
Tri-channel	NP	●	Unigrip™
	RP	●	

**Table 2 – Healing Screw – Compatible Implant Platforms and Screwdrivers**

Slim Healing Abutment for	Available platforms	Color coding	Screwdriver
Tri-channel	NP	●	Unigrip™
	RP	●	

**Table 3 – Slim Healing Abutments – Compatible Implant Platforms and Screwdrivers**

Healing Abutment Anatomical PEEK for	Available platforms	Color coding	Screwdriver
Conical Connection (CC)	WP	None	Unigrip™

**Table 4 – Healing Abutments Anatomical PEEK – Compatible Implant Platforms and Screwdrivers**

Healing abutments which feature tri-oval conical connection and Healing Abutments Anatomical PEEK are co-packed with a clinical screw. Refer to Nobel Biocare Instructions for Use (IFU) IFU1057 for information regarding clinical screws. This IFU is available for download at [ifu.nobelbiocare.com](http://ifu.nobelbiocare.com).

## Intended Use / Intended Purpose

### Healing Abutments / Healing Screw

Intended to be temporarily connected to an endosseous dental implant or implant abutment, to support healing of the surrounding soft tissue.

## Indications

Healing abutments and Healing screws are indicated for use with endosseous dental implants or implant abutments in the maxilla or mandible for supporting single tooth to full arch denture procedures.

Healing abutments Nobel Biocare N1™ TCC are indicated for use for up to 180 days.

Healing Abutments Bridge Conical Connection are additionally indicated to prevent growth of bone on the implant platform to support placement of an Impression Coping Bridge.

Healing Abutments Anatomical PEEK Conical Connection are adjustable abutments which are indicated for use for up to 180 days with endosseous dental implants in the maxilla or mandible, for supporting single tooth to full arch denture procedures.

## Contraindications

TI is contraindicated to use healing abutments in:

- Patients who are medically unfit for an oral surgical procedure.
- Patients in whom adequate sizes, numbers or desirable positions of implants are not reachable to achieve safe support of functional or eventually parafunctional loads.
- Patients who are allergic or hypersensitive to titanium alloy Ti-6Al-4V (titanium, aluminum, vanadium), DLC (Diamond Like Carbon) or PEEK (Polyetheretherketone).

## Cautions

### General

Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for successful implant treatment.

It is strongly recommended that Nobel Biocare surgical instruments and prosthetic components are used only with Nobel Biocare implants, as combining components that are not dimensioned for correct mating can lead to mechanical and/or instrumental failure, damage to tissue or unsatisfactory esthetic results.

### Before Surgery

Special attention must be given to patients who have local or systemic factors that could interfere with the healing process of either bone or soft tissue or the osseointegration process (e.g. cigarette smoking, poor oral hygiene, uncontrolled diabetes, oro-facial radiotherapy, steroid therapy, infections in the neighboring bone). Special caution is advised in patients who receive bisphosphonate therapy.

In general, implant placement and prosthetic design must accommodate individual patient conditions. In case of bruxism, other parafunctional habits or unfavorable jaw relationships, reappraisal of the treatment option may be considered.

The device has not been evaluated in pediatric/adolescent patients and is not recommended for use in children. Routine treatment is not recommended for pediatric patients until the end of the jawbone growth phase has been properly documented.

Pre-operative hard tissue or soft tissue deficits may yield a compromised esthetic result or unfavorable implant angulations.

All components, instruments and tooling used during the clinical and laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

## At Surgery

Because of the small size of the devices, care must be taken that they are not swallowed or aspirated by the patient. It is appropriate to use specific supporting tools to prevent aspiration of loose parts (e.g. gauze, a dental dam or a throat shield).

## After Surgery

The Screwdriver Manual Nobel Biocare N1™ Base and Screwdriver Machine Nobel Biocare N1™ Base are used to tighten and loosen the clinical screw which fastens the Nobel Biocare N1™ Base to the dental implant.

## Intended Users and Patient Groups

Healing abutments are to be used by dental health care professionals.

Healing abutments are to be used in patients subject to dental implant treatment.

## Clinical Benefits and Undesirable Side Effects

### Clinical Benefits Associated with Healing Abutments

Healing abutments are a component of treatment with a dental implant system and/or dental crowns and bridges. As a clinical benefit of treatment, patients can expect to have their missing teeth replaced and/or crowns restored.

### Undesirable Side Effects Associated with Healing Abutments

The placement of these devices is part of an invasive treatment which may be associated with typical side effects such as inflammation, infection, bleeding, hematoma, pain, and swelling. During abutment placement or removal the pharyngeal reflex (gag reflex) may be triggered in patients with a sensitive gag reflex.

Healing abutments are part of a multi-component system that replaces teeth and as a result, the implant recipient may experience side effects similar to those associated with teeth, such as retained cement, calculus, mucositis, ulcers soft tissue hyperplasia, soft and/or hard tissue recessions. Some patients may experience discoloration in the mucosal area such as graying.

## Notice regarding serious incidents

For a patient / user / third party in the European Union and in countries with an identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. The contact information for the manufacturer of this device to report a serious incident is as follows:

Nobel Biocare AB  
[www.nobelbiocare.com/complaint-form](http://www.nobelbiocare.com/complaint-form)

## Handling Procedure

### Handling Procedure for Healing Abutments, Healing Screw and Slim Healing Abutments

1. Select appropriate abutment and check occlusal clearance.
2. Connect the abutment to implant and hand-tighten using dedicated screwdriver. See Table 1 or Table 2 for compatibility. Refer to Nobel Biocare IFU1085 for information regarding the screwdriver. It is recommended to verify the final abutment seating using radiographic imaging.

**Caution** Never exceed recommended tightening torque for the screw. Overtightening of abutment screw may lead to a screw fracture.

3. If removal of the abutment is needed, untighten it using dedicated screwdriver.
4. For abutments featuring tri-oval conical connection, if the removal is not possible, use the Abutment Retrieval tool. Refer to Nobel Biocare IFU1096 for information regarding the Abutment Retrieval Tool.

### Handling Procedure for Healing Abutments Anatomical PEEK

1. Select appropriate abutment and check occlusal clearance.
2. If necessary, adjust the abutment height using a rotary instrument (e.g. carbide or acrylic bur). The tissue facing axial contours of the abutment may be modified to achieve the desired shape. If axial modification is done, polishing with silicone points or similar methods is recommended.
3. Connect to implant and hand-tighten using Unigrip™ screwdriver. Refer to Nobel Biocare IFU1085 for information regarding the screwdriver. It is recommended to verify the final abutment seating using radiographic imaging.
4. If removal of the abutment is needed, untighten the screw using dedicated screwdriver.
5. If the abutment cannot be removed, use the Abutment Retrieval tool. Refer to Nobel Biocare IFU1096 for information regarding the Abutment Retrieval Tool.

## Materials

- Healing Abutments, Healing Screw and Slim Healing Abutments: Titanium alloy 90% Ti, 6% Al, 4% V according to ASTM F136 and ISO 5832-3.
- Healing Abutments Anatomical PEEK: PEEK (Polyetheretherketone).
- Clinical screws for Healing Abutments Nobel Biocare N1™ TCC and for Healing Abutments Anatomical PEEK: Titanium alloy 90% Ti, 6% Al, 4% V according to ASTM F136 and ISO 5832-3 and DLC (Diamond Like Carbon) coating.

# Sterility and Reusability Information

Healing abutments have been sterilized using irradiation and are intended for single use. Do not use after the labeled expiration date.

**Warning** Do not use device if the packaging has been damaged or previously opened as the device sterility and/or integrity may be compromised.

**Caution** Healing abutments are single use products and must not be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause local or systemic infection.

# Magnetic Resonance (MR) Safety Information

## MR Safety Information for Single Tooth Restoration

MRI Safety Information



Non-clinical testing has demonstrated that Healing abutments are MR conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions mentioned here below. Failure to follow these conditions may result in injury to the patient.

Nominal value(s) of Static Magnetic Field [T]	1.5-Tesla (1.5 T)	3-Tesla (3 T)
Maximum Spatial Field Gradient [T/m and gauss/cm]	Maximum spatial field gradient of 58.9 T/m (5,890 G/cm).	
RF Excitation	Circularly Polarized (CP).	
RF Transmit Coil Type	Whole body transmit coil.	
Maximum Whole-Body SAR [W/kg]	Inferior to the neck: 2.0 W/kg Superior to the neck: 0.5 W/kg	Inferior to the xyphoid: 2.0 W/kg Between the xyphoid and neck: 1.0 W/kg Superior to the neck: 0.5 W/kg
Limits on Scan Duration	Under the scan conditions defined above, the dental implant systems are expected to produce a maximum temperature rise less than 6.0°C after 15 minutes of continuous scanning.	
MR Image Artifact	In non-clinical testing, the image artifact caused by the dental implant systems extend radially approximately 3.0 cm from the devices or device assemblies when imaged in a 3 T MRI system.	

## MR Safety Information for Single Tooth Configurations Zygoma Implants (applicable only during Zygoma Implant healing phase)

MRI Safety Information



Non-clinical testing has demonstrated the Brånemark System® Zygoma Healing Abutment, Healing Abutment Bmk Syst RP (when used with Zygoma) is MR conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions mentioned here below. Failure to follow these conditions may result in injury to the patient.

Nominal value(s) of Static Magnetic Field [T]	1.5-Tesla (1.5 T)	3-Tesla (3 T)
Maximum Spatial Field Gradient [T/m and gauss/cm]	Maximum spatial field gradient of 58.9 T/m (5,890 G/cm).	
RF Excitation	Circularly Polarized (CP).	
RF Transmit Coil Type	Whole body transmit coil.	
Maximum Whole-Body SAR [W/kg]	Inferior to the shoulder: 2.0 W/kg Superior to the shoulder: 0.2 W/kg	Inferior to the xyphoid: 2.0 W/kg Superior to the xyphoid: 0.2 W/kg

Limits on Scan Duration	Under the scan conditions defined above, the dental implant systems are expected to produce a maximum temperature rise less than 6.0°C after 15 minutes of continuous scanning.
MR Image Artifact	In non-clinical testing, the image artifact caused by the dental implant systems extend radially approximately 2.4 cm from the devices or device assemblies when imaged in a 3 T MRI system.
Caution	Configurations with more than 2 Zygoma implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of configurations with more than 2 Zygoma implants in the MR environment is unknown. Scanning a patient who has this configuration may result in patient injury.

Implant placement with intention to restore at prosthetic level with PIBs or IBOs (multiple tooth restorations):

Please consult IFU for NobelProcera® Implant Bridge Titanium and Zirconia, NobelProcera® Crown and Bridge, NobelProcera® HT ML FCZ, and NobelProcera® Implant Bar Overdenture for use as part of a bridge configuration.

# Performance Requirements and Limitations

To achieve the desired performance, healing abutments must only be used with the products described in this Instructions for Use and/or in the Instructions for Use for other compatible Nobel Biocare products, and in accordance with the Intended Use for each product. To confirm the compatibility of products which are intended to be used in conjunction with healing abutments, check the color coding, dimensions, lengths, connection type and/or any direct marking as applicable on the products or product labeling.

# Facilities and Training

It is strongly recommended that new and experienced users of Nobel Biocare products always go through special training before using a new product for the first time. Nobel Biocare offers a wide range of courses for various levels of knowledge and experience. For more information please visit [www.nobelbiocare.com](http://www.nobelbiocare.com).

# Storage, Handling and Transportation




The device must be stored and transported in dry conditions in the original packaging at room temperature and not exposed to direct sunlight. Incorrect storage and transportation may influence device characteristics leading to failure.

# Disposal

Safely discard potentially contaminated or no longer usable medical devices as healthcare (clinical) waste in accordance with local healthcare guidelines, country and government legislation or policy.

Separation, re-cycling or disposal of packaging material shall follow local country and government legislation on packaging and packaging waste, where applicable.

# Manufacturer and Distributor Information

	<p>Nobel Biocare AB PO Box 5190, 402 26 Västra Hamngatan 1 411 17 Göteborg Sweden <a href="http://www.nobelbiocare.com">www.nobelbiocare.com</a></p>
<p>UK Responsible Person</p> <div style="border: 1px solid black; padding: 2px; display: inline-block;">UK RP</div>	<p>Nobel Biocare UK Ltd 4 Longwalk Road Stockley Park Uxbridge UB11 1FE United Kingdom</p>
<p>Distributed in Australia by</p>	<p>Nobel Biocare Australia Pty Ltd Level 4, 7 Eden Park Drive Macquarie Park, NSW 2113 Australia Phone: +61 1800 804 597</p>
<p>Distributed in New Zealand by</p>	<p>Nobel Biocare New Zealand Ltd 33 Spartan Road Takanini, Auckland, 2105 New Zealand Phone: +64 0800 441 657</p>
<p>CE Mark for Class IIb Devices</p>	
<p>UKCA Mark for Class IIb Devices</p>	

**Note** Regarding Canadian Device Licensure not all products described in the IFU may have a device licence according to Canadian law.

## Basic UDI-DI Information

The following table lists the Basic UDI-DI information for the devices described in this IFU.



















































Product	Basic UDI-DI Number
Healing Abutments Conical Connection 3.0/NP/RP/WP	73327470000001236T
Healing Abutments Bridge Conical Connection NP/RP/WP	
Healing Abutments Nobel Biocare N1™ TCC NP/RP	
Healing Abutments NobelReplace® NP/RP/WP/6.0	
Healing Abutments Brånemark System® NP/RP/WP	
Brånemark System® Zygoma Healing Abutments	
Slim Healing Abutments Conical Connection 3.0/NP/RP	
Healing Abutments Anatomical PEEK Conical Connection WP	

EN All rights reserved.

Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Product images in this folder are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product.

# Symbols Glossary

The following symbols may be present on the device labeling or in information accompanying the device. Refer to the device labeling or accompanying information for the applicable symbols.

							
Authorized Representative in the European Community/ European Union	UK Responsible Person	Authorised Representative in Switzerland	Sterilized using Ethylene Oxide	Sterilized using irradiation	Sterilized using steam or dry heat		
							
Batch code	Catalogue number	Unique Device Identifier	Serial number	Medical device	Magnetic resonance safe		
							
Caution	Magnetic resonance conditional	Non-sterile	Contains hazardous substances	Contains or presence of DEHP phthalate	Contains or presence of natural rubber latex	Contains or presence of phthalate	Contains biological material of animal origin
						 <a href="https://symbol.glossary.nobelbiocare.com">symbol.glossary.nobelbiocare.com</a> <a href="https://ifu.nobelbiocare.com">ifu.nobelbiocare.com</a>	
CE mark	CE mark with Notified Body number	UKCA mark	UKCA mark with Approved Body number	Consult instructions for use	For prescription use only	Link to Online Symbols Glossary and IFU Portal	
							
Date of manufacture	Manufacturer	Use-by date	Upper limit of temperature	Temperature limit	Do not resterilize	Do not re-use	Non-pyrogenic
							
Date	Tooth number	Patient number	Patient identification	Health care centre or doctor	Patient information website	EU Importer	Swiss Importer
							
Double sterile barrier system	Single sterile barrier system	Single sterile barrier system with protective packaging inside	Single sterile barrier system with protective packaging outside	Do not use if package is damaged and consult instructions for use	Keep away from sunlight	Keep dry	