creos™

mucogain

ENGLISH	4	SRPSKI	106
DEUTSCH	10	ČEŠTINA	112
FRANÇAIS	16	ROMÂNĂ	118
ESPAÑOL	22	TÜRKÇE	124
PORTUGUÊS	28		
ITALIANO	34		
SVENSKA	40		
DANSKA	46		
SUOMI	52		
NORSK	58		
РУССКИЙ	64		
LIETUVIŲ	70		
NEDERLANDS	76		
MAGYAR	82		
POLSKI	88		
ΕΛΛΗΝΙΚΑ	94		
HRVATSKI	100		

## CREOS™ MUCOGAIN COLLAGEN MATRIX

### DESCRIPTION

creos™ mucogain is a resorbable matrix intended to provide support during coverage procedures of localized gingival recessions and for local soft tissue augmentation. creos™ mucogain is designed to provide an off-the-shelf alternative to autogenous soft tissue grafts.

The open, interconnecting porous structure of creos™ mucogain provides a matrix for the migration of proliferating cells and vascular structures. As the healing process advances, creos™ mucogain is degraded whilst a new soft tissue matrix is regenerated inside the creos™ mucogain structure.

creos™ mucogain is composed of highly purified porcine collagen and elastin fibers. Starting materials for the production of creos™ mucogain are harvested in EU certified facilities under strict veterinary controls from animals declared

fit for human consumption. creos™ mucogain is produced from collagen using controlled and validated manufacturing processes and is not chemically crosslinked. creos™ mucogain is double packaged and sterilized by ethylene oxide gas treatment.

### PROPERTIES

creos™ mucogain has a unique open porous structure. This structure was designed to guarantee a good regeneration as a result of migration of cells and blood vessels as well as the initial nutrient supply. creos™ mucogain is biodegradable and will be resorbed during the healing process whilst in parallel a new soft tissue matrix is generated.

The excellent biocompatibility of creos™ mucogain is due to the utilization of effective purification and processing methods during manufacturing and has been demonstrated both in vitro as well as in vivo. creos™ mucogain can be fixed in place using sutures or pins. creos™ mucogain is sterile and non-pyrogenic. Experimental

studies have shown that creos™ mucogain can be considered safe with regard to potential viral transmission.

Since creos™ mucogain has good hemostatic properties, the risk of haematoma occurring during treatment is diminished. creos™ mucogain is hydrophilic and retains its structural integrity also in the wet state.

### INDICATIONS FOR USE

creos™ mucogain is intended to be used for soft tissue augmentation indications in the oral cavity around teeth or implants:

- Guided tissue regeneration (GTR) procedures in recession defects for root coverage.
- Localized gingival augmentation to increase keratinized tissue around teeth and implants.

### LIMITATIONS OF USE

Contraindications:

- creos™ mucogain is not indicated for exposed healing.

Adverse Reactions:

- Adverse reactions during the use of collagen based products during dental soft tissue regeneration procedures have not been observed. However, allergic reactions towards porcine collagen may not be totally excluded.

Precautions:

- In patients with constrained healing capacity (e.g. due to metabolic disorders, cancer treatments, heavy smoking) creos™ mucogain should be used with special caution and enhanced postoperative care.
- creos™ mucogain should not be used in cases where active infection or inflammation exists. In such cases, the infection or inflammation should be adequately treated first.

- creos™ mucogain is not recommended for use in patients known to be sensitive to materials of porcine origin.
- creos™ mucogain should not be used to increase keratinized tissue around teeth and implants if no remaining keratinized tissue is present in the recipient area. Therefore, the surgical preparation and treatment should be performed in a way that a minimal amount of keratinized tissue is preserved in the recipient area.
- The defect is exposed by means of (a) properly prepared flap(s) and the usual surgical procedures in creating a properly prepared site are undertaken.
- creos™ mucogain can be shaped to the required size with surgical scissors or surgical scalpels. The trimming can be performed in either dry or wet state.
- creos™ mucogain is applied over the prepared site and stabilized using resorbable sutures or the flap. Excessive pressure should be avoided as it may compress the matrix.
- Complete penetration of creos™ mucogain by blood and exudates allows close adaptation and adhesion of the matrix to the underlying surface and facilitates the formation of a blood clot within creos™ mucogain.
- Due to the tensile strength of this matrix, fixation or suturing is possible. Fixation or suturing of the matrix may be indicated, depending upon the nature of the particular defect and to avoid dislocation of the matrix. The flap should be sutured over the matrix without tension.

#### INSTRUCTIONS FOR USE

The principles of sterile handling, appropriate patient medication and the general surgical techniques must be followed.

Root coverage procedures should be preceded by a conditioning with cause-related anti-infective therapy to eradicate any remaining bacterial infection and a thorough training of the patient in oral hygiene. During treatment the root surface should be debrided and planed carefully.

#### POST-OPERATIVE CARE

The usual postoperative care and medication should be given. Further prosthetic treatment should only be carried out after an initial healing period to ensure complete soft tissue regeneration. Previous studies have indicated that initial soft tissue healing, i.e., wound closure and the resolution of normal inflammation, occurs within 4 – 8 weeks. But clinical judgment, taking into account the healing of both soft and hard tissues in the treated patient, should be applied before prosthetic treatment. In case of bacterial contamination rinsing with appropriate bactericidal solutions is recommended.

In the event that early matrix removal is necessary, the tissues adjacent to the matrix should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual matrix. Following careful separation of the surrounding tissue, the remaining portion of the matrix can be excised and the area debrided to remove any inflamed or infected tissue.

#### WARNINGS AND OTHER INFORMATION

Do not use creos™ mucogain when the sterile packaging is damaged.

creos™ mucogain cannot be resterilized and is for single use only. A resterilization will lead to a deterioration of product properties and may result in unwanted biological reactions due to an unsterile product or to a treatment failure due to an early loss of the matrix stability.

#### STORAGE AND HANDLING

creos™ mucogain should be stored at room temperature ( $\leq 25^{\circ}\text{C}$  /  $\leq 77^{\circ}\text{F}$ ) and in a dry place. creos™ mucogain should be handled using sterile gloves and/or with sterile instruments. Do not use after the use by date indicated on the packaging.

The manufacturer is not liable for an unintended or incorrect use of creos™ mucogain.

**PRESENTATION FORM**

One unit of creos™ mucogain is provided per double blister package.

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Single use

Consult instructions for use

Caution. Consult post-operative care and limitations of use

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Matricel GmbH
Nobel Biocare AB, Box 5190
Matricel GmbH, Kaiserstrasse 100
Distributor: Nobel Biocare AB, Box 5190, 402 26,