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IFU1007. Brånemark System® Mk III TiUnite®, NobelSpeedy®, Twist Drill, Twist Step Drill, Counterbore, Screw Tap

Nome prodotto	Vite di copertura	Driver per impianto	Viti	Abutment			Transfer per impronta
				Guarigione	Provvisorio	Permanente	
Brånemark System® Mk III TiU NP	Cover Screw Brånemark System® NP	Implant Driver NP Brånemark System Implant Driver Wrench Adapter NP	Abutment Screw Multi-unit Angled Prosthetic Screw Multi-unit	Healing Abutment NP	Immediate Temporary Abutment NP	Procera Esthetic Abutment NP	Impression Coping Closed Tray NP
					Plastic Coping Immediate NP	Esthetic Abutment NP	Impression Coping Open Tray NP
					Temporary Abutment NP	15° Esthetic Abutment NP	Impression Coping Open Tray Multi-unit NP
					QuickTemp™ Abutment Conical NP	Snappy™ Abutment NP	
					Universal Base NP	Universal Base NP	
					GoldAdapt™ Engaging NP	GoldAdapt™ Engaging NP	Impression Coping Closed Tray Multi-unit NP
					Plastic Coping QuickTemp™ Abutment Conical NP	GoldAdapt™ Non- Engaging NP	
					Temporary Abutment Engaging NP	Multi-Unit Abutment NP	Impression Coping Bar Closed Tray Multi- unit NP
					Temporary Abutment Non-Engaging NP	17° Multi-Unit Abutment NP	Temporary Coping Multi-unit Titanium NP
					Temporary Abutment Plastic Engaging NP		Temporary Coping Multi-unit Plastic NP
					Temporary Abutment Plastic Non-Engaging NP		
					Brånemark System® Mk III TiU RP	Cover Screw Brånemark System® RP	Implant Driver RP Brånemark System® Implant Driver Wrench Adapter RP
Plastic Coping Immediate Temporary Abutment RP	Esthetic Abutment RP	Impression Coping Open Tray RP					
QuickTemp™ Abutment Conical RP	15° Esthetic Abutment RP	Impression Coping Bar Brånemark System® RP					
Universal Base RP	Snappy™ Abutment RP						
Universal Base Burn-out Coping RP	Universal Base RP						
Plastic Coping QuickTemp™ Abutment Conical RP	GoldAdapt™ Engaging RP	Impression Coping Open Tray Multi-unit RP					
Temporary Abutment Engaging RP	GoldAdapt™ Non- Engaging RP	Impression Coping Closed Tray Multi- unit RP					
Temporary Abutment Non-Engaging RP	Multi-Unit Abutment RP	Impression Coping Bar Closed Tray Multi- unit RP					
Temporary Abutment Plastic Engaging RP	17° Multi-Unit Abutment RP	Temporary Coping Multi-unit Titanium RP					
Temporary Abutment Plastic Non-Engaging RP	30° Multi-Unit Abutment RP	Temporary Coping Multi-unit Plastic RP					
	30° Multi-Unit Abutment Non-Engaging RP						
	Locator Abutment RP						
	Ball Abutment Titanium RP						
Brånemark System® Mk III TiU WP	Cover Screw Brånemark System® WP	Implant Driver WP Brånemark System Implant Driver Wrench Adapter WP	Prosthetic Screw Multi-unit	Healing Abutment WP	Immediate Temporary Abutment WP	Procera Esthetic Abutment WP	Impression Coping Closed Tray WP
					Plastic Coping Immediate Temporary Abutment WP	Esthetic Abutment WP	Impression Coping Open Tray WP
					QuickTemp™ Abutment Conical WP	15° Esthetic Abutment WP	Impression Coping Open Tray Multi-unit WP
					Universal Base WP	Snappy™ Abutment WP	
					Universal Base Burn-out Coping WP	Universal Base WP	
					Plastic Coping QuickTemp™ Abutment Conical WP	GoldAdapt™ Engaging WP	Impression Coping Closed Tray Multi-unit WP
					Temporary Abutment Engaging WP	GoldAdapt™ Non- Engaging WP	Impression Coping Bar Closed Tray Multi-unit WP
					Temporary Abutment Non-Engaging WP	Multi-Unit Abutment WP	Temporary Coping Multi-unit Titanium WP
					Temporary Abutment Plastic Engaging WP	Locator Abutment WP	Temporary Coping Multi-unit Plastic WP
					Temporary Abutment Plastic Non-Engaging WP		

Nome prodotto	Vite di copertura	Driver per impianto	Viti	Abutment			Transfer per impronta
				Guarigione	Provvisorio	Permanente	
NobelSpeedy® Shorty NP	Cover Screw Brånemark System® NP	Implant Driver NP Brånemark System® Implant Driver Wrench Adapter NP	Abutment Screw Multi-unit Angled Prosthetic Screw Multi-unit	Healing Abutment NP	Immediate Temporary Abutment NP	Procera Esthetic Abutment NP	Impression Coping Closed Tray NP
NobelSpeedy® Groovy NP					Plastic Coping Immediate NP Temporary Abutment NP QuickTemp™ Abutment Conical NP Plastic Coping QuickTemp™ Abutment Conical NP Temporary Abutment Engaging NP Temporary Abutment Non-Engaging NP Temporary Abutment Plastic Engaging NP Temporary Abutment Plastic Non-Engaging NP	Esthetic Abutment NP 15° Esthetic Abutment NP Snappy™ Abutment NP Universal Base NP Universal Base Burn-out Coping NP GoldAdapt™ Engaging NP GoldAdapt™ Non-Engaging NP Multi-Unit Abutment NP 17° Multi-Unit Abutment NP Locator Abutment NP Gold Coping Bar Multi-unit (Prosthetic screw included) NP	Impression Coping Open Tray NP Impression Coping Open Tray Multi-unit NP Impression Coping Closed Tray Multi-unit NP Impression Coping Bar Closed Tray Multi-unit NP Temporary Coping Multi-unit Titanium NP Temporary Coping Multi-unit Plastic NP
NobelSpeedy® Shorty RP	Cover Screw Brånemark System® RP	Implant Driver RP Brånemark System Implant Driver Wrench Adapter RP		Healing Abutment RP	Immediate Temporary Abutment RP	Procera Esthetic Abutment RP	Impression Coping Closed Tray RP
NobelSpeedy® Groovy RP					Plastic Coping Immediate Temporary Abutment RP QuickTemp™ Abutment Conical RP Plastic Coping QuickTemp™ Abutment Conical RP Temporary Abutment Engaging RP Temporary Abutment Non-Engaging RP Temporary Abutment Plastic Engaging RP Temporary Abutment Plastic Non-Engaging RP	Esthetic Abutment RP 15° Esthetic Abutment RP Snappy™ Abutment RP Universal Base RP Universal Base Burn-out Coping RP GoldAdapt™ Engaging RP GoldAdapt™ Non-Engaging RP Multi-Unit Abutment RP 17° Multi-Unit Abutment RP 30° Multi-Unit Abutment RP 30° Multi-Unit Abutment Non-Engaging RP Locator Abutment RP Ball Abutment Titanium RP Gold Abutment Bar Implant Level (vite clinica inclusa) RP Gold Coping Bar Multi-unit (vite protesica inclusa) RP	Impression Coping Open Tray RP Impression Coping Bar Brånemark System® RP Impression Coping Open Tray Multi-unit RP Impression Coping Closed Tray Multi-unit RP Impression Coping Bar Closed Tray Multi-unit RP Temporary Coping Multi-unit Titanium RP Temporary Coping Multi-unit Plastic RP
NobelSpeedy® Shorty WP	Cover Screw Brånemark System® WP	Implant Driver WP Brånemark System Implant Driver Wrench Adapter WP	Prosthetic Screw Multi-unit	Healing Abutment WP	Immediate Temporary Abutment WP	Procera Esthetic Abutment WP	Impression Coping Closed Tray WP
NobelSpeedy® Groovy WP					Plastic Coping Immediate Temporary Abutment WP QuickTemp™ Abutment Conical WP Plastic Coping QuickTemp™ Abutment Conical WP Temporary Abutment Engaging WP Temporary Abutment Non-Engaging WP Temporary Abutment Plastic Engaging WP Temporary Abutment Plastic Non-Engaging WP	Esthetic Abutment WP 15° Esthetic Abutment WP Snappy™ Abutment WP Universal Base WP Universal Base Burn-out Coping WP GoldAdapt™ Engaging WP GoldAdapt™ Non-Engaging WP Multi-Unit Abutment WP Locator Abutment WP	Impression Coping Open Tray WP Impression Coping Open Tray Multi-unit WP Impression Coping Closed Tray Multi-unit WP Impression Coping Bar Closed Tray Multi-unit WP Temporary Coping Multi-unit Titanium WP Temporary Coping Multi-unit Plastic WP

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IFU1008. NobelReplace® Select™ TC

Nome prodotto	Vite di copertura	Driver per impianto	Viti cliniche	Abutment			Transfer per impronta	Altro
				Guarigione	Provvisorio	Permanente		
Replace Select™ TC NP	Healing Screw RplSel TC NP Cover Screw NobRpl NP	Implant Driver NobRpl NP	Screw Multi-unit Angled Abutm NobRpl NP Screw Ceramic Abutment NobRpl NP Abutment Screw NobRpl NP	Healing Abutment NobRpl NP	Temporary Abutment Engaging NobRpl NP Temporary Abutment Non-Eng NobRpl NP Temporary Abutment Plastic Eng NobRpl NP Temp Abutm Plastic Non-Eng NobRpl NP Immediate Temporary Abutment NobRpl NP	Snappy™ Abutments NobRpl NP Esthetic Abutment NobRpl NP 15° Esthetic Abutment NobRpl NP Universal Base Tri-Channel NP GoldAdapt™ Engaging NobRpl NP GoldAdapt™ Non-Engaging NobRpl NP Multi-unit Abutment NobRpl NP 17° Multi-unit Abutment NobRpl NP Ball Abutment Titanium NobRpl NP NobelProcera® Ab Ti NobelReplace® NP w-u NobelProcera® Ab Ti NobelReplace® NP Procera® Abutm Zr NobRpl NP	Impr Coping Open Tray NobRpl NP Ø 4,5 mm Impr Coping Closed Tray NobRpl NP Impr Cop Cl Tray Low Prof NRpl NP Impr Coping Open Tray NobRpl NP	PIB IBO Impl Retr Instr Hex & Tri-Ch NP/RP 22 mm Bone Mill with Guide NobRpl NP Ø 4,6 Implant Rescue Collar Tri-Channel Ø 3,5
Replace Select™ TC RP	Healing Screw RplSel TC RP Cover Screw NobRpl RP	Implant Driver NobRpl RP	Screw Multi-unit Angled Abutm NobRpl RP Screw Ceramic Abutment NobRpl RP/WP/6.0 Abutment Screw NobRpl RP/WP/6.0	Healing Abutment NobRpl RP	Temporary Abutment Engaging NobRpl RP Temporary Abutment Non-Eng NobRpl RP Temporary Abutment Plastic Eng NobRpl RP Temp Abutm Plastic Non-Eng NobRpl RP Immediate Temporary Abutment NobRpl RP	Snappy™ Abutments NobRpl RP Esthetic Abutment NobRpl RP 15° Esthetic Abutment NobRpl RP Universal Base Tri-Channel RP GoldAdapt™ Engaging NobRpl RP GoldAdapt™ Non-Engaging NobRpl RP Multi-unit Abutment NobRpl RP 30° Multi-unit Abutment NobRpl RP 30° Multi-u Abutm Non-Eng NobRpl RP 17° Multi-unit Abutment NobRpl RP Ball Abutment Titanium NobRpl NP NobelProcera® Ab Ti NobelReplace® RP w-u NobelProcera® Ab Ti NobelReplace® RP Procera® Abutm Zr NobRpl RP	Impr Coping Closed Tray NobRpl RP Impr Cop Cl Tray Low Prof NRpl RP Impr Coping Open Tray NobRpl RP	PIB IBO Impl Retr Instr Hex & Tri-Ch NP/RP 31 mm Bone Mill with Guide NobRpl RP Ø 5,3 Implant Rescue Collar Tri-Channel Ø 4,3

IFU1010. NobelReplace® Conical Connection, NobelReplace® Conical Connection PMC, NobelReplace® Conical Connection TiUltra™

Nome prodotto	Vite di copertura	Driver per impianto	Viti cliniche	Abutment			Transfer per impronta	Altro
				Guarigione	Provvisorio	Permanente		
NobelReplace® Conical Connection NP	Cover Screw CC NP	Implant driver CC NP	Clinical Screw CC NP	Healing abutment CC NP	Slim Temporary Abutment CC NP	Universal Base CC NP	Open Tray CC NP	Guided Implant Mount NobRpl CC NP Bone Mill Guide CC NP Guide Pin CC NP NobelProcera®: NobelProcera® Abutment Zirconia CC NP NobelProcera® Abutment Ti CC NP NobelProcera® FCZ Implant Crown ASC CC NP Implant Bridge Restorations (Ti, Zr) Fixed e Fixed Removable Bars (Ti) Metal Adapter CC NP
NobelReplace® Conical Connection PMC NP			Omnigrip™ Clinical Screw CC NP	Healing abutment bridge NP	Temporary Snap Abutment CC NP	Esthetic Abutment CC NP	Closed Tray CC NP	
NobelReplace® Conical Connection TiUltra™ NP			On1™ Clinical Screw NP Screw MU Angled Abutment CC NP	Slim healing abutment CC NP	Temporary Abutment Non Engaging CC NP	Procera Esthetic Abutment CC NP Multi-Unit Abutment Plus CC NP (0°/17°/30°) On1™ Base NP	Bridge Open Tray CC NP	
NobelReplace® Conical Connection RP	Cover Screw CC RP	Implant driver CC RP	Clinical Screw CC RP/WP	Healing abutment CC RP	Slim Temporary Abutment CC RP	Universal Base CC NP	Open Tray CC RP	Guided Implant Mount NobRpl CC RP Bone Mill Guide CC RP Guide Pin CC RP NobelProcera®: NobelProcera® ASC Abutment Zirconia CC RP NobelProcera® Abutment Ti CC RP NobelProcera® FCZ Implant Crown ASC CC RP Implant Bridge Restorations (Ti, Zr) Fixed e Fixed Removable Bars (Ti) Metal Adapter CC RP
NobelReplace® Conical Connection PMC RP			Omnigrip™ Clinical Screw CC RP/WP	Healing abutment bridge RP	Temporary Snap Abutment CC RP	Esthetic Abutment CC RP	Closed Tray CC RP	
NobelReplace® Conical Connection TiUltra™ RP			On1™ Clinical Screw RP Screw MU Angled Abutment CC RP	Slim healing abutment CC RP	Temporary Abutment Non Engaging CC RP	Procera Esthetic Abutment CC RP Multi-Unit Abutment Plus CC RP (0°/17°/30°) On1™ Base RP	Bridge Open Tray CC RP	

IFU1012. Replace Select™ Tapered TiUnite®, Replace Select™ Tapered PMC

Nome prodotto	Vite di copertura	Driver per impianto	Viti cliniche	Abutment			Transfer per impronta	Altro
				Guarigione	Provvisorio	Permanente		
Replace Select™ Tapered NP	Cover Screw NobRpl NP Healing Screw RplSel TC NP	Implant Driver NobRpl NP	Abutment Screw NobelReplace® NP	Healing Abutment NobRpl NP	Temporary Abutment Engaging NobRpl NP	Snappy Abutment NobRpl NP Esthetic Abutment NobRpl NP	Impr Coping Open Tray NobRpl NP	PIB IBO Bone Mill with Guide NobRpl NP Bone Mill Guide NobelReplace® NP Impl Retr Instr Hex & Tri-Ch NP/RP Implant Rescue Collar Tri-Channel Ø 3,5
			Screw Multi-unit Angled Abutment NobRpl NP Screw Ceramic Abutment NobRpl NP	Temporary Abutment Non-Eng NobRpl NP Immediate Temporary Abutment NobRpl NP Plastic/Temp Coping Esthetic Abutment NP Plastic/Temp Coping 15° Esthetic Ab NP	15° Esthetic Abutment NobRpl NP Universal Base Tri-Channel NP Universal Base Burn-out Coping NP/RP GoldAdapt™ Engaging NobRpl NP GoldAdapt™ Non-Engaging NobRpl NP Multi-unit Abutment NobRpl NP 17° Multi-unit Abutment NobRpl NP Locator Abut for 3.5 RPL Sel Impl. Ball Abutment Titanium NobRpl NP NobelProcera® Ab Ti NobelReplace® NP Procera® Abutment Zr NobRpl NP			
Replace Select™ Tapered RP	Cover Screw NobRpl RP Healing Screw RplSel TC RP	Implant Driver NobRpl RP	Abutment Screw NobRpl RP/WP/6.0	Healing Abutment NobRpl RP	Temporary Abutment Engaging NobRpl RP	Snappy Abutment NobRpl RP Esthetic Abutment NobRpl RP	Impl Retr Instr Hex & Tri-Ch NP/RP	PIB IBO Bone Mill Guide NobelReplace® RP Bone Mill with Guide NobRpl RP Impl Retr Instr Hex & Tri-Ch NP/RP Implant Rescue Collar Tri-Channel Ø 4,3
			Screw Multi-unit Angled Abutment NobRpl RP Screw Ceramic Abutment NobRpl RP/WP/6.0	Temporary Abutment Non-Eng NobRpl RP Immediate Temporary Abutment NobRpl RP	15° Esthetic Abutment NobRpl RP Universal Base Tri-Channel RP Universal Base Burn-out Coping NP/RP GoldAdapt™ Engaging NobRpl RP GoldAdapt™ Non-Engaging NobRpl RP Multi-unit Abutment NobRpl RP 17° Multi-unit Abutment NobRpl RP 30° Multi-unit Abutment NobRpl RP Locator Abut for 4.3 RPL Sel Impl. Ball Abutment Titanium NobRpl RP 3 mm NobelProcera® Ab Ti NobelReplace® RP Procera® Abutment Zr NobRpl RP			

Nome prodotto	Vite di copertura	Driver per impianto	Viti cliniche	Abutment			Transfer per impronta	Altro
				Guarigione	Provvisorio	Permanente		
Replace Select™ Tapered WP	Cover Screw NobRpl WP	Implant Driver NobRpl WP	Screw Ceramic Abutment NobRpl RP/ WP/6.0	Healing Abutment NobRpl WP	Temporary Abutment Engaging NobRpl WP Temporary Abutment Non-Eng NobRpl WP Temp Abutment Plastic Non-Eng NobRpl WP Immediate Temporary Abutment NobRpl WP Plastic/Temp Cop Esthetic Ab NobRpl WP Pl/Temp Coping 15° Esthetic Ab NobRpl WP	Snappy Abutment NobRpl WP Esthetic Abutment NobRpl WP 15° Esthetic Abutment NobRpl WP Universal Base Burn-out Coping WP Universal Base Tri-Channel WP GoldAdapt™ Non-Engaging NobRpl WP GoldAdapt™ Engaging NobRpl WP Multi-unit Abutment NobRpl WP Locator Abut for RPL Sel Impl. NobelProcera® Ab Ti NobelReplace® WP Procera® Abutment Zr NobRpl WP	Impr Coping Open Tray NobRpl WP	PIB IBO Bone Mill Guide NobelReplace® WP Bone Mill with Guide NobRpl WP Ø 6,5 Impl Retr Instr CC RP & Tri-Ch WP Impl Retr Instr CC RP& Tri-Ch WP& TCC RP
Replace Select™ Tapered 6.0	Cover Screw NobRpl 6.0 mm	Implant Driver NobRpl 6.0	Screw Ceramic Abutment NobRpl RP/ WP/6.0	Healing Abutment NobRpl 6.0	Temporary Abut Plastic Eng NobRpl 6.0 Temp Abutment Plastic Non-Eng NobRpl 6.0 Temporary Abutment Engaging NobRpl 6.0 Temporary Abutment Non-Eng NobRpl 6.0 Immediate Temporary Abutment NobRpl 6.0	Snappy Abutment NobRpl 6.0 Esthetic Abutment NobRpl 6.0 15° Esthetic Abutment NobRpl 6.0 GoldAdapt™ Non-Engaging NobRpl 6.0 GoldAdapt™ Engaging NobRpl 6.0 NobelProcera® Ab Ti NobelReplace® 6.0 Procera® Abutment Zr NobRpl 6.0	Impr Coping Open Tray NobRpl 6.0	PIB IBO Bone Mill Guide NobRpl 6.0 Bone Mill with Guide NobRpl 6.0 Ø 7 Impl Retr Instr Tri-Channel 6.0 22 mm

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Vite di copertura	Tipo di connessione	Piattaforma implantare	Sistema implantare	Cacciavite
Cover Screw Brånemark System® NP	Esagonale esterna	NP	Brånemark System® NobelSpeedy® Groovy	Cover Screw Driver Brånemark System Hexagon (IFU1085)
Cover Screw Brånemark System® RP		RP	Brånemark System® NobelSpeedy® Groovy NobelZygoma™ 0°	
Cover Screw Brånemark System® WP		WP	Brånemark System® NobelSpeedy® Groovy	
Cover Screw Conical Connection 3.0	Connessione conica	3.0	NobelActive®	Unigrip™ (IFU1085)
Cover Screw Conical Connection NP		NP	NobelActive® NobelParallel™ CC NobelReplace® CC	
Cover Screw Conical Connection RP		RP	NobelActive® NobelParallel™ CC NobelReplace® CC Sistemi Trefoil™	
Cover Screw Conical Connection WP		WP	NobelActive® NobelParallel™ CC NobelReplace® CC	
Cover Screw Nobel Biocare N1™ TCC NP	Connessione conica triovale	NP	Nobel Biocare N1™ TiUltra™ TCC	Omnigrip™ Mini (IFU1085)
Cover Screw Nobel Biocare N1™ TCC RP		RP	Nobel Biocare N1™ TiUltra™ TCC	
Cover Screw NobelReplace® NP	Interna trilobata	NP	NobelReplace® Tapered Groovy, ReplaceSelect™	Unigrip™ (IFU1085)
Cover Screw NobelReplace® RP		RP		
Cover Screw NobelReplace® WP		WP		
Cover Screw NobelReplace® 6.0 mm		6.0		

IFU1018. Esthetic Abutment Conical Connection/Brånemark System®/ NobelReplace®, Esthetic Abutment Nobel Biocare N1™ TCC

Esthetic Abutment	Tipo di connessione	Vite clinica	Codifica a colori	Cacciavite	Cappetta in plastica/provvisoria
Esthetic Abutment Conical Connection 3.0	Connessione conica	Clinical Screw Conical Connection 3.0	Nessuno	Unigrip™	-
Esthetic Abutment Conical Connection NP	Connessione conica	Clinical Screw Conical Connection NP	●	Unigrip™	-
Esthetic Abutment Conical Connection RP	Connessione conica	Clinical Screw Conical Connection RP/ WP	●	Unigrip™	-
Esthetic Abutment Conical Connection WP	Connessione conica	Clinical Screw Conical Connection RP/ WP	●	Unigrip™	-
Esthetic Abutment Bmk System NP	Esagonale esterna	Abutment Screw Brånemark System® NP	Nessuno	Unigrip™	Plastic/ Temporary Coping Esthetic Abutment
Esthetic Abutment Bmk System RP	Esagonale esterna	Abutment Screw Brånemark System® RP	Nessuno	Unigrip™	Plastic/ Temporary Coping Esthetic Abutment
Esthetic Abutment Bmk System WP	Esagonale esterna	Abutment Screw Brånemark System® WP	Nessuno	Unigrip™	Plastic/ Temporary Coping Esthetic Abutment
Esthetic Abutment NobelReplace® NP	Trilobata	Abutment Screw NobelReplace® NP	●	Unigrip™	Plastic/ Temporary Coping Esthetic Abutment
Esthetic Abutment NobelReplace® RP	Trilobata	Abutment Screw NobelReplace® RP/WP/6.0	●	Unigrip™	Plastic/ Temporary Coping Esthetic Abutment

Esthetic Abutment	Tipo di connessione	Vite clinica	Codifica a colori	Cacciavite	Cappetta in plastica/provisoria
Esthetic Abutment NobelReplace® WP	Trilobata	Abutment Screwtrub NobelReplace® RP/WP/6.0	●	Unigrip™	Plastic/Temporary Coping Esthetic Abutment
Esthetic Abutment NobelReplace® 6.0	Trilobata	Abutment Screwtrub NobelReplace® RP/WP/6.0	●	Unigrip™	Plastic/Temporary Coping Esthetic Abutment
Esthetic Abutment Nobel Biocare N1™ TCC NP	Connessione conica triovale	Clinical Screw Nobel Biocare N1 TCC NP	● (vite)	Omnigrip™ Mini	-
Esthetic Abutment Nobel Biocare N1™ TCC RP	Connessione conica triovale	Clinical Screw Nobel Biocare N1 TCC RP	● (vite)	Omnigrip™ Mini	-

IFU1023. Universal Base e Universal Abutment

Universal Base/ Abutment per	Piattaforme disponibili	Altezze di margine disponibili	Torque di serraggio	Cacciavite
Connessione conica (CC)	NP RP WP	1,5 mm 3,0 mm	35 Ncm	Unigrip™
Connessione conica trilobata (TCC)	NP RP	1,5 mm 3,0 mm	20 Ncm	Omnigrip™ mini
Trilobata	NP RP WP	1,5 mm 3,0 mm	35 Ncm	Unigrip™
Esagonale esterna	NP RP WP	1,5 mm 3,0 mm	35 Ncm	Unigrip™

IFU1040. Abutment Release Pin CC

Abutment	Abutment Release Pin CC 3.0	Abutment Release Pin CC NP	Abutment Release Pin CC RP/WP
Temporary Abutment Engaging Conical Connection	✓	n/d	n/d
Esthetic Abutment Conical Connection	✓	n/d	n/d
GoldAdapt™ Engaging Conical Connection	✓	n/d	n/d
NobelProcera® Abutment Zirconia Conical Connection	n/d	✓	✓

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IFU1048. Frese Tapered Drill e Dense Bone Drill, maschiatori Screw Tap

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Vite clinica/ dell'abutment/ protesica	Abutment/Manufatto	Cacciavite
Abutment Screw Brånemark System®	Universal Base External Hex	Unigrip™
	Esthetic Abutment Brånemark System®	
	Temporary Abutment Brånemark System®	
	NobelProcera® Abutment Ti Brånemark System®	
	NobelProcera® Implant Bridge Ti Brånemark System®	
	NobelProcera® Implant Bar Overdenture Brånemark System®	
	Snappy™ Abutment Brånemark System® GoldAdapt™ Brånemark System®	
	Gold Abutment Brånemark System®	
Screw Ceramic Abutment Brånemark System®	NobelProcera® Abutment Zr Brånemark System®	Unigrip™
	NobelProcera® HT ML FCZ Implant Bridge Brånemark System®	
	NobelProcera® Zirconia Implant Bridge Brånemark System®	
	Procera® Esthetic Abutment Brånemark System®	
Screw Multi-unit Angled Abutment CC	17°/30° Multi-unit Abutment CC	Unigrip™
Vite clinica Multi-unit Abutment Nobel Biocare N1™ TCC	Straight Multi-unit Abutment Xael™ TCC	Multi-unit
Screw Multi-unit Angled Abutment NobelReplace®	17°/30° Multi-unit Abutment NobelReplace®	Unigrip™
Screw Multi-unit Angled Abutment Brånemark System®	17°/30° Multi-unit Abutment Brånemark System®	Unigrip™

Vite clinica/ dell'abutment/ protesica	Abutment/Manufatto	Cacciavite
Vite protesica Multi-unit Abutment	Temporary Coping Multi-unit Abutment NobelProcera® Implant Bridge Ti Multi-unit	Unigrip™
	NobelProcera® Implant Bridge Zr Multi-unit	
	NobelProcera® Implant Bar Overdenture Multi-unit	
	NobelProcera® HT ML FCZ Implant Bridge Multi-unit	
	Gold Coping Multi-unit	
Prosthetic Screw Multi-unit Abutment Omnigrip™ Mini	NobelProcera® Zirconia Implant Bridge Multi-unit	Omnigrip™ Mini
NobelZygoma 0° Angled Multi-unit Abutment Screw	45°/60° Multi-unit Abutment External Hex	Unigrip™
Brånemark System® Zygoma Abutment Screw	NobelProcera® Implant Bridge Zygoma	Unigrip™
Brånemark System® Zygoma Angled Multi-unit Abutment Screw	17° Multi-unit Abutment Brånemark System® Zygoma	Unigrip™
Omnigrip™ Clinical Screw Titanium (CC & Tri-Channel)	Nobel Procera® Titanium Abutment ASC (CC & Tri-Channel)	Omnigrip™

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IFU1074. Concetto On1™

Nome prodotto	Impianto	Abutment	Vite	Driver
On1™ Base CC NP 1.75 mm	Impianti con connessione conica interna: - NobelActive® - NobelActive® TiUltra™ - NobelReplace® CC - NobelReplace® CC TiUltra™ - NobelParallel™ CC - NobelParallel™ CC TiUltra™	On1™ Temporary Abutment Engaging	On1™ Clinical Screw 1.75 mm NP	On1™ Screwdriver
On1™ Base CC NP 2.5 mm		On1™ Temporary Abutment Non-Engaging	On1™ Clinical Screw 2.5 mm NP	
On1™ Base CC RP 1.75 mm		On1™ Universal Abutment	On1™ Clinical Screw 1.75 mm RP	
On1™ Base CC RP 2.5 mm		On1™ Universal Abutment Non-Engaging	On1™ Clinical Screw 2.5 mm RP	
On1™ Base CC WP 1.75 mm		On1™ Esthetic Abutment	On1™ Clinical Screw 1.75 mm WP	
On1™ Base CC WP 2.5 mm			On1™ Clinical Screw 2.5 mm WP	
On1™ Universal Abutment 0.3 mm NP		N/D	On1™ Base	
On1™ Universal Abutment 1.25 mm NP	On1™ Base Xael		On1™ Prosthetic Lab Screw	
On1™ Universal Abutment 0.3 mm RP				
On1™ Universal Abutment 1.25 mm RP				
On1™ Universal Abutment 0.3 mm WP				
On1™ Universal Abutment 1.25 mm WP				
On1™ Esthetic Abutment Titanium 0.3 mm NP	N/D	On1™ Base	On1™ Prosthetic Screw	Unigrip™ Screwdriver
On1™ Esthetic Abutment Titanium 0.3 mm RP		On1™ Base Xael	On1™ Prosthetic Lab Screw	
On1™ Esthetic Abutment Titanium 0.3 mm WP				
On1™ Base Xael™ CC NP 1.75 mm	Impianti con connessione conica interna: - NobelActive® - NobelActive® TiUltra™ - NobelReplace® CC - NobelReplace® CC TiUltra™ - NobelParallel™ CC - NobelParallel™ CC TiUltra™	On1™ Temporary Abutment Engaging	On1™ Clinical Screw 1.75 mm NP	On1™ Screwdriver
On1™ Base Xael™ CC NP 2.5 mm		On1™ Temporary Abutment Non-Engaging	On1™ Clinical Screw 2.5 mm NP	
On1™ Base Xael™ CC RP 1.75 mm		On1™ Universal Abutment	On1™ Clinical Screw 1.75 mm RP	
On1™ Base Xael™ CC RP 2.5 mm		On1™ Universal Abutment Non-Engaging	On1™ Clinical Screw 2.5 mm RP	
On1™ Base Xael™ CC WP 1.75 mm		On1™ Esthetic Abutment	On1™ Clinical Screw 1.75 mm WP	
On1™ Base Xael™ CC WP 2.5 mm			On1™ Clinical Screw 2.5 mm WP	

Nome prodotto	Impianto	Abutment	Vite	Driver
On1™ Universal Abutment Non-Eng NP 0.3 mm	N/D	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw On1™ Prosthetic Lab Screw	Unigrip™ Screwdriver
On1™ Universal Abutment Non-Eng NP 1.25 mm				
On1™ Universal Abutment Non-Eng RP 0.3 mm				
On1™ Universal Abutment Non-Eng RP 1.25 mm				
On1™ Universal Abutment Non-Eng WP 0.3 mm				
On1™ Universal Abutment Non-Eng WP 1.25 mm				
On1™ Healing Cap NP 1.5 mm	N/D	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw	Unigrip™ Screwdriver
On1™ Healing Cap NP 2.5 mm				
On1™ Healing Cap RP 1.5 mm				
On1™ Healing Cap RP 2.5 mm				
On1™ Healing Cap WP 1.5 mm				
On1™ Healing Cap WP 2.5 mm				
On1™ IOS Healing Cap 6 mm NP	N/D	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw	Unigrip™ Screwdriver
On1™ IOS Healing Cap 4.5 mm NP				
On1™ IOS Healing Cap 6 mm RP				
On1™ IOS Healing Cap 4.5 mm RP				
On1™ IOS Healing Cap 5 mm WP				
On1™ IOS Healing Cap 4 mm WP				
On1™ Temporary Abutment Engaging NP	N/D	On1™ Base On1™ Base Xeal	On1™ Prosthetic Screw On1™ Prosthetic Lab Screw	Unigrip™ Screwdriver
On1™ Temporary Abutment Non-Engaging NP				
On1™ Temporary Abutment Engaging RP				
On1™ Temporary Abutment Non-Engaging RP				
On1™ Temporary Abutment Engaging WP				
On1™ Temporary Abutment Non-Engaging WP				
On1™ Impression Coping Closed Tray NP	N/D	On1™ Base On1™ Base Xeal	N/D	Unigrip™ Screwdriver
On1™ Impression Coping Closed Tray RP				
On1™ Impression Coping Closed Tray WP				
On1™ Impression Coping Open Tray NP				
On1™ Impression Coping Open Tray RP				
On1™ Impression Coping Open Tray WP				
On1™ Impression Coping Open Tray N-Eng NP	N/D	On1™ Base On1™ Base Xeal	N/D	Unigrip™ Screwdriver
On1™ Impression Coping Open Tray N-Eng RP				
On1™ Impression Coping Open Tray N-Eng WP				
On1™ Clinical Screw 1.75 mm NP	Impianti con connessione conica interna: - NobelActive® - NobelActive® TiUltra™ - NobelReplace® CC - NobelReplace® CC TiUltra™ - NobelParallel™ CC - NobelParallel™ CC TiUltra™	On1™ Temporary Abutment Engaging	On1™ Prosthetic Screw	On1™ Screwdriver
On1™ Clinical Screw 2.5 mm NP		On1™ Temporary Abutment Non-Engaging	On1™ Prosthetic Lab Screw	
On1™ Clinical Screw 1.75 mm RP		On1™ Universal Abutment	On1™ Impression Coping Screw	
On1™ Clinical Screw 2.5 mm RP		On1™ Universal Abutment Non-Engaging		
On1™ Clinical Screw 1.75 mm WP		On1™ Esthetic Abutment		
On1™ Clinical Screw 2.5 mm WP				
On1™ Prosthetic Screw	N/D		On1™ Clinical Screw	Unigrip™ Screwdriver
On1™ Screwdriver Machine 21 mm	N/D	On1™ Base On1™ Base Xeal	On1™ Clinical Screw	N/D
On1™ Screwdriver Manual 25 mm	N/D	On1™ Base On1™ Base Xeal	On1™ Clinical Screw	N/D

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Connection to Handpiece

Connection to Handpiece	Dispositivi compatibili	Linea di impianti
Connection to Handpiece	Guided Implant Mount Brånemark System® NP, RP 3.75, RP 4.0, WP Zygoma Implant Mount	Brånemark System® NobelSpeedy® Zygoma RP
NobelReplace® Connection to Handpiece	Guided Implant Mount NobelActive® NP, RP 4.3, RP 5.0, WP 5.5 Guided Implant Mount NobelParallel™ NP 3.75, RP 4.3, RP 5.0, WP 5.5 Guided Implant Mount NobelReplace® NP 3.5, RP 4.3, WP 5.0, 6.0 Guided Implant Mount NobelReplace® CC NP 3.5, RP 4.3, RP 5.0 Implant Mount Ø 5 mm	NobelActive® NobelParallel™ NobelReplace® Tapered NobelReplace® CC Sistema Trefoil™

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Direction Indicator Ø 2/Ø 2.4-2.8 mm e Direction Indicator Tapered NP/RP/WP/6.0

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Handle for Implant Rescue Collar e Drill Guide

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Implant Driver Wrench Adapter Brånemark System® NP/RP/WP

Impianto	Istruzioni per l'uso	Torque Wrench
Brånemark System®	IFU1007	Brånemark System® Manual Torque Wrench
NobelSpeedy®	IFU1007	

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
Nome prodotto	Dispositivi compatibili		
	Drill Stop/Drill Stop Kit Box	Frese	Cacciaviti
Drill Stop Ø 2.8	Guided Drill Stop Kit Box Drill Stop Kit Box	Guided Twist Drill 2.8x(10+)7-13 mm Guided Twist Drill 2.8x(10+)7-18 mm Twist Step Drill 2.4/2.8 10-18 mm Twist Step Drill 2.4/2.8 7-15 mm Twist Step Drill 2.4/2.8 7-10 mm Guided Twist Step Drill Ø 2.4/2.8 (10+) 7-13 mm Guided Twist Step Drill Ø 2.4/2.8 (10+) 7-18 mm Twist Step Drill 2.4/2.8x18-25 mm, LS Twist Step Drill 2.4/2.8x18-25 mm Guided Drill Guide NP per Ø 2,8 mm Guided Drill Guide RP per Ø 2,8 mm Guided Drill Guide 6.0/WP per Ø 2,8 mm	Screwdriver Machine Unigrip™ Screwdriver Manual Unigrip™
Drill Stop Ø 3	Guided Drill Stop Kit Box Drill Stop Kit Box	Guided Twist Drill 3x(10+)7-13 mm Guided Twist Drill 3x(10+)7-18 mm Twist Drill 3x10-18 mm Twist Drill 3x7-15 mm Twist Drill 3x7-10 mm Guided Drill Guide NP per Ø 3 mm Guided Drill Guide RP per Ø 3 mm Guided Drill Guide 6.0/WP per Ø 3 mm	Screwdriver Machine Unigrip™ Screwdriver Manual Unigrip™
Drill Stop Ø 3.2	Guided Drill Stop Kit Box Drill Stop Kit Box	Twist Step Drill 2.8/3.2 7-10 mm Twist Drill 3.2x18-25 mm Twist Step Drill 2.8/3.2 10-18 mm Twist Step Drill 2.8/3.2 7-15 mm Guided Twist Drill 3.2x(10+)7-13 mm Twist Drill 3.2x7-10 mm Guided Twist Drill 3.2x(10+)7-18 mm Twist Drill 3.2x10-18 mm Guided Twist Step Drill Ø 2.8/3.2 (10+) 7-18 mm Guided Twist Step Drill Ø 2.8/3.2 (10+) 7-13 mm Twist Drill 3.2x7-15 mm Guided Drill Guide NP per Ø 3.2 mm Guided Drill Guide RP per Ø 3.2 mm	Screwdriver Machine Unigrip™ Screwdriver Manual Unigrip™
Drill Stop Ø 3.4	Guided Drill Stop Kit Box Drill Stop Kit Box	Guided Twist Drill 3.4x(10+)7-13 mm Guided Twist Drill 3.4x(10+)7-18 mm Twist Drill 3.4x10-18 mm Twist Drill 3.4x7-15 mm Twist Drill 3.4x7-10 mm Twist Drill 3.4x18-25 mm Guided Drill Guide RP per Ø 3,4 mm	Screwdriver Machine Unigrip™ Screwdriver Manual Unigrip™
Drill Stop Ø 3.6	Guided Drill Stop Kit Box Drill Stop Kit Box	Twist Step Drill 3.2/3.6 7-15 mm Twist Step Drill 3.2/3.6 7-10 mm Guided Twist Step Drill Ø 3.2/3.6 (10+) 7-18 mm Guided Twist Step Drill Ø 3.2/3.6 (10+) 7-13 mm Twist Step Drill 3.2/3.6 10-18 mm Guided Drill Guide RP per Ø 3.6 mm Guided Drill Guide WP/6.0 per Ø 3.6 mm	Screwdriver Machine Unigrip™ Screwdriver Manual Unigrip™
Drill Stop Ø 3.8	Guided Drill Stop Kit Box Drill Stop Kit Box	Guided Twist Drill 3.8x(10+)7-13 mm Guided Twist Drill 3.8x(10+)7-18 mm Guided Drill Guide 6.0/WP per Ø 3.8 mm	Screwdriver Machine Unigrip™ Screwdriver Manual Unigrip™
Drill Stop Ø 4.2	Guided Drill Stop Kit Box Drill Stop Kit Box	Guided Twist Step Drill 3.8/4.2 (10+) 7-18 mm Twist Step Drill 3.8/4.2 10-18 mm Twist Step Drill 3.8/4.2 7-15 mm Twist Step Drill 3.8/4.2 7-10 mm Guided Twist Drill 4.2x(10+)7-13 mm Guided Twist Drill 4.2x(10+)7-18 mm Guided Twist Step Drill 3.8/4.2 (10+) 7-13 mm Guided Drill Guide RP per Ø 4.2 mm Guided Drill Guide 6.0/WP per Ø 4.2 mm	Screwdriver Machine Unigrip™ Screwdriver Manual Unigrip™

Nome prodotto	Dispositivi compatibili		
	Drill Stop/Drill Stop Kit Box	Frese	Cacciaviti
Drill Stop Kit Box	Drill Stop Ø 2 Drill Stop Ø 2.8 Drill Stop Ø 3 Drill Stop Ø 3.2 Drill Stop Ø 3.4 Drill Stop Ø 3.6 Drill Stop Ø 3.8 Drill Stop Ø 4.2	Twist Drill con punta 2.0x18-25 mm Twist Drill con punta 2x7-15 mm Twist Drill 3.2x7-10 mm Twist Drill 1.5x7-15 mm Twist Drill con punta 2x10-18 mm Twist Drill con punta 2x18-25 mm, LS Twist Step Drill 2.4/2.8 10-18 mm Twist Step Drill 2.4/2.8 7-15 mm Twist Step Drill 2.4/2.8x18-25 mm, LS Twist Step Drill 2.4/2.8x18-25 mm Twist Step Drill 2.4/2.8x18-25 mm Twist Drill 3x10-18 mm Twist Drill 3x7-15 mm Twist Drill 3x7-10 mm Twist Step Drill 2.8/3.2 7-10 mm Twist Drill 3.2x18-25 mm Twist Step Drill 2.8/3.2 10-18 mm Twist Step Drill 2.8/3.2 7-15 mm Twist Drill 3.2x7-10 mm Twist Drill 3.2x10-18 mm Twist Drill 3.2x7-15 mm Twist Drill 3.4x10-18 mm Twist Drill 3.4x7-15 mm Twist Drill 3.4x7-10 mm Twist Drill 3.4x18-25 mm Twist Step Drill 3.2/3.6 7-15 mm Twist Step Drill 3.2/3.6 7-10 mm Twist Step Drill 3.2/3.6 10-18 mm Twist Step Drill 3.8/4.2 10-18 mm Twist Step Drill 3.8/4.2 7-15 mm Twist Step Drill 3.8/4.2 7-10 mm	N/D
Guided Drill Stop Kit Box	Drill Stop Ø 2 Drill Stop Ø 2.8 Drill Stop Ø 3 Drill Stop Ø 3.2 Drill Stop Ø 3.4 Drill Stop Ø 3.6 Drill Stop Ø 3.8 Drill Stop Ø 4.2	Guided Twist Drill 2x(10+)7-13 mm Guided Twist Drill 2x(10+)7-18 mm Guided Twist Drill 2.8x(10+)7-13 mm Guided Twist Step Drill Ø 2.4/2.8 (10+) 7-13 mm Guided Twist Step Drill Ø 2.4/2.8 (10+) 7-18 mm Guided Drill Guide NP per Ø 2.8 mm Guided Drill Guide RP per Ø 2.8 mm Guided Drill Guide 6.0/WP per Ø 2.8 mm Guided Twist Drill 3x(10+)7-13 mm Guided Twist Drill 3x(10+)7-18 mm Guided Twist Drill 3.2x(10+)7-13 mm Guided Twist Drill 3.2x(10+)7-18 mm Guided Twist Step Drill Ø 2.8/3.2 (10+) 7-18 mm Guided Twist Step Drill Ø 2.8/3.2 (10+) 7-13 mm Guided Twist Drill 3.4x(10+)7-13 mm Guided Twist Drill 3.4x(10+)7-18 mm Guided Twist Step Drill Ø 3.2/3.6 (10+) 7-18 mm Guided Twist Step Drill Ø 3.2/3.6 (10+) 7-13 mm Guided Twist Drill 3.8x(10+)7-13 mm Guided Twist Drill 3.8x(10+)7-18 mm Guided Twist Step Drill 3.8/4.2 (10+) 7-18 mm Guided Twist Drill 4.2x(10+)7-13 mm Guided Twist Step Drill 3.8/4.2 (10+) 7-13 mm	N/D

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N. articolo	Nome prodotto	Vite di copertura	Implant Mount e relativa vite	Transfer per impronta e relativo perno guida	Abutment di guarigione	Abutment definitivi diritti e relative viti	Abutment definitivi angolati e relativa vite	Ponti su impianti e Implant Bar Overdenture e relative viti	Zygoma Handle, Connection to Handpiece Altro
28862	Zygoma Implant RP 30 mm	Zygoma Implant Cover Screw : 28989 Brånemark System Zygoma Implant Cover Screw : 32424	Zygoma Implant Mount Zygoma Fixture Mount Screw M2specx6,6 mm (il montaimpianto è preassemblato all'impianto)	Brånemark System® Zygoma Impression Coping Open Tray Ø 4 mm Guide Pin Zygoma 20 mm RP	Brånemark System® Zygoma Healing Abutment Ø 4x3 mm Brånemark System® Zygoma Healing Abutment Ø 4x5 mm	Zygoma Abutment Multi-unit RP 3 mm Zygoma Abutment Multi-unit RP 5 mm Brånemark System® Zygoma Multi-unit Abutment 3 mm Brånemark System® Zygoma Multi-unit Abutment 5 mm	Zygoma 17° Abutment Multi-unit RP 2 mm Zygoma 17° Abutment Multi-unit RP 3 mm Brånemark System® Zygoma 17° Multi-unit Abutment 2 mm Brånemark System® Zygoma 17° Multi-unit Abutment 3 mm Brånemark System® Zygoma Angled Multi-unit Abutment Screw	Nobel Procera Implant Bar Overdenture Titanium Procera Implant Bridge Ti livello dell'impianto NobelProcera® Zr Implant Bridge (livello dell'impianto) Brånemark System® Zygoma Abutment Screw	Zygoma Handle Connection to Handpiece Unigrip Screw Driver (alla vite di montaggio dell'impianto)
28863	Zygoma Implant RP 35 mm								
28864	Zygoma Implant RP 40 mm								
28865	Zygoma Implant RP 42,5 mm								
28866	Zygoma Implant RP 45 mm								
28867	Zygoma Implant RP 47,5 mm								
28868	Zygoma Implant RP 50 mm								
28869	Zygoma Implant RP 52,5 mm								

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Strumento	Diametro dell'impianto compatibile
Trepine Drill 3.2/4.0 mm	3,0 mm
Trepine Drill 3.8/4.6 mm	3,3 mm/3,5 mm/3,75 mm
Trepine Drill 4.4/5.2 mm	4,0 mm/4,3 mm
Trepine Drill 5.2/6.2 mm	5,0 mm
Trepine Drill 5.6/6.6 mm	5,5 mm
Trepine Drill 6.2/7.0 mm	6,0 mm

IFU1098. Manual Torque Wrench Surgical e Prosthetic

Chiave	Adattatore per chiave	Strumenti
NobelActive® Manual Torque Wrench Surgical	NOBELREPLACE® Manual Torque Wrench Adapter Surgical	Implant Driver CC 3.0 28 mm (IFU1090)
		Implant Driver CC 3.0 per Slim Abutment Implant Driver CC NP (IFU1090) Implant Driver CC NP per Slim Abutment (IFU1090) Implant Driver CC RP (IFU1090)
NobelReplace® Manual Torque Wrench Surgical	NOBELREPLACE® Manual Torque Wrench Adapter Surgical	Implant Driver CC RP per Slim Abutment (IFU1090)
		Implant Driver CC WP (IFU1090) Strumenti per la rimozione dell'impianto (IFU1097) Abutment Screw Retrieval Instrument (IFU1043)
Brånemark System® Manual Torque Wrench Surgical	Brånemark System® Manual Torque Wrench Surgical	Implant Driver Brånemark Syst NP (IFU1090)
		Implant Driver Brånemark Syst RP (IFU1090) Implant Driver Brånemark Syst WP (IFU1090)
Manual Torque Wrench Surgical Nobel Biocare N1™	NobelReplace® Manual Torque Wrench Adapter Surgical	Implant Driver Nobel Biocare N1™ TCC NP (IFU1087)
		Implant Driver Nobel Biocare N1™ TCC RP (IFU1087) Strumenti per la rimozione dell'impianto (IFU1097) Abutment Screw Retrieval Instrument (IFU1043)
Manual Torque Wrench Prosthetic Nobel Biocare N1™	Manual Torque Wrench Adapter Prosthetic	Screwdriver Machine Unigrip (IFU1085)
		Screwdriver Machine Multi-Unit (IFU1085) Omnigrip™ Screwdriver Machine (IFU1085) On1™ Base Screwdriver Machine (IFU1074) Screwdriver Machine Multi-unit Brånemark System (IFU1085) Screwdriver Machine Ball Abutment (IFU1085)
Manual Torque Wrench Prosthetic Nobel Biocare N1™	Manual Torque Wrench Adapter Prosthetic	Screwdriver Machine Omnigrip™ Mini (IFU1085)
		Screwdriver Machine Multi-Unit (IFU1085)

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Abutment originale	Vite sostitutiva	Torque di serraggio	Cacciavite sostitutivo
Standard Abutment RP o EsthetiCone Abutment	Prosthetic Screw Internal Hexagon	10 Ncm	Screwdriver Hexagon 27 mm Machine Screwdriver Hex Long
	Prosthetic Screw Slot		Screwdriver Medium 37 mm Screwdriver Slot Short 27 mm
	Prosthetic Screw Conical		Machine Screwdriver Slot Short Machine Screwdriver Slot Long
Brånemark 3.0 NP Abutment per 3.0 Brånemark Implant	Converter Screw Titanium Unigrip™ fit Ø 3.0	15 Ncm	Unigrip™ Screwdriver*

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Abutment originale	Cacciaviti sostitutivi
Ball Abutment	Screwdriver Machine Ball Abutment Screwdriver Manual Ball Abutment

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IFU2009. Guided Pilot Drill Sleeve, Guided Sleeve e Guided Anchor Pin Sleeve

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Nome prodotto	Impianto	Manicotto	Cacciavite	Altro
Guided Drill Guide 6.0/WP Ø 6 per Ø 5 mm	N/D	Guided Sleeve 6.0/WP	N/D	Handle for Guided Drill Guide
Guided Drill Guide 6.0/WP per Ø 2,8 mm				
Guided Drill Guide 6.0/WP per Ø 2 mm				
Guided Drill Guide 6.0/WP per Ø 3,8 mm				
Guided Drill Guide 6.0/WP per Ø 3 mm				
Guided Drill Guide 6.0/WP per Ø 4,2 mm				
Guided Drill Guide NP per Ø 2,8 mm		Guided Sleeve NP		
Guided Drill Guide NP per Ø 2 mm				
Guided Drill Guide NP per Ø 3,2 mm				
Guided Drill Guide NP per Ø 3 mm				
Guided Drill Guide RP per Ø 3,2 mm		Guided Sleeve RP		
Guided Drill Guide RP per Ø 2,8 mm				
Guided Drill Guide RP per Ø 2 mm				
Guided Drill Guide RP per Ø 3,4 mm				
Guided Drill Guide RP per Ø 3,6 mm				
Guided Drill Guide RP per Ø 3 mm				
Guided Drill Guide RP per Ø 4,2 mm				
Guided Drill Guide Tapered RP per NP				
Guided Drill Guide Tapered 6.0/WP per NP		Guided Sleeve 6.0/WP		
Guided Drill Guide Tapered 6.0/WP per RP				
Guided Drill Guide WP/6.0 per Ø 3,6 mm				
Guided Drill Guide WP/6.0 per Ø 4,6 mm				
Guided Surgery Implant Mount NobelActive® NP	NobelActive®/TiUltra™ NP	Guided Sleeve NP	Surgical Driver	Manual Torque Wrench Adapter Surgical Connection to Handpiece
Guided Surgery Implant Mount NobelActive® RP 4.3	NobelActive®/TiUltra™ RP 4.3	Guided Sleeve RP	Screwdriver Manual/Machine Unigrip™	
Guided Surgery Implant Mount NobelActive® RP 5.0	NobelActive®/TiUltra™ RP 5.0	Guided Sleeve 6.0/WP		
Guided Surgery Implant Mount NobelActive® WP 5.5	NobelActive®/TiUltra™ WP 5.5			

Nome prodotto	Impianto	Manicotto	Cacciavite	Altro
Guided Implant Mount Bmk Syst NP	Brånemark System® NP	Guided Sleeve NP	Surgical Driver	Brånemark System® Manual Torque Wrench Adapter Surgical Connection to handpiece
Guided Implant Mount Bmk Syst RP	Brånemark System® RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	
Guided Implant Mount Bmk Syst WP	Brånemark System® WP	Guided Sleeve 6.0/ WP		
Guided Implant Mount NobRpl NP	NobelReplace® NP	Guided Sleeve NP	Surgical Driver	NobelReplace® Manual Torque Wrench Adapter Surgical Connection to handpiece
Guided Implant Mount NobRpl RP	NobelReplace® RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	
Guided Implant Mount NobRpl WP	NobelReplace® WP	Guided Sleeve 6.0/ WP	Guided Sleeve RP	
Guided Implant Mount NobRpl 6.0	NobelReplace® 6.0	Guided Sleeve 6.0/ WP	Guided Sleeve 6.0/ WP	
Guided Implant Mount NobelReplace® Conical Connection NP 3.5	NobelReplace® CC NP	Guided Sleeve NP	Surgical Driver	NobelReplace® Manual Torque Wrench Adapter Surgical Connection to handpiece
Guided Implant Mount NobelReplace® Conical Connection RP 4.3	NobelReplace® CC RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	
Guided Implant Mount NobelReplace® Conical Connection RP 5.0	NobelReplace® CC RP 5.0	Guided Sleeve 6.0/ WP		
Guided Implant Mount NobelParallel™ CC NP	NobelParallel™ CC NP	Guided Sleeve NP	Surgical Driver	NobelReplace® Manual Torque Wrench Adapter Surgical Connection to handpiece
Guided Implant Mount NobelParallel™ CC RP 4.3	NobelParallel™ CC RP	Guided Sleeve RP	Screwdriver Manual/ Machine Unigrip™	
Guided Implant Mount NobelParallel™ CC RP 5.0	NobelParallel™ CC RP 5.0	Guided Sleeve 6.0/ WP		
Guided Implant Mount NobelParallel™ CC WP 5.5	NobelParallel™ CC WP	Guided Sleeve 6.0/ WP		

Nome prodotto	Impianto	Manicotto	Cacciavite	Altro
Guided Template Abutm w Screw NobelActive® NP	NobelActive®/ NobelActive® TiUltra™ NP	Guided Sleeve NP	Screwdriver Manual/ Machine Unigrip™	N/D
Guided Template Abutm w Screw NobelActive® RP 4.3	NobelActive®/ NobelActive® TiUltra™ RP 4.3	Guided Sleeve RP		
Guided Template Abutm w Screw NobelActive® RP 5.0	NobelActive®/ NobelActive® TiUltra™ RP 5.0	Guided Sleeve 6.0/ WP		
Guided Template Abutm w Screw NobRpl 6.0	NobelReplace® CC/ NobelReplace® CC TiUltra™ 6.0 e Replace Select™ Tapered TiUnite® 6.0 implants			
Guided Template Abutm w Screw NobRpl NP	NobelReplace® CC/ NobelReplace® CC TiUltra™ NP e Replace Select™ Tapered TiUnite® NP implants	Guided Sleeve NP		
Guided Template Abutm w Screw NobRpl RP	NobelReplace® CC/ NobelReplace® CC TiUltra™ RP e Replace Select™ Tapered TiUnite® RP implants	Guided Sleeve RP		
Guided Template Abutm w Screw NobRpl WP	NobelReplace® CC/ NobelReplace® CC TiUltra™ WP e Replace Select™ Tapered TiUnite® WP implants	Guided Sleeve 6.0/ WP		
Guided Template Abutment CC NP	NobelActive®/ NobelActive® TiUltra™ NP, NobelParallel™ CC/ NobelParallel™ CC TiUltra™ NP, NobelReplace® CC/ NobelReplace® CC TiUltra™, NP NobelReplace® CC PMC NP implants	Guided Sleeve NP		
Guided Template Abutment CC NP	NobelReplace® CC/ NobelReplace® CC TiUltra™ & CC PMC NP implants			
Guided Template Abutment CC RP 4.3	NobelReplace® CC/ NobelReplace® CC TiUltra™ & CC PMC RP 4.3 implants	Guided Sleeve RP		
Guided Template Abutment CC RP 5.0	NobelReplace® CC/ NobelReplace® CC TiUltra™ & CC PMC RP 5.0 implants	Guided Sleeve 6.0/ WP		
Guided Template Abutment CC WP 5.5	NobelReplace® CC/ NobelReplace® CC TiUltra™ & CC PMC WP 5.5 implants			
Guided Template Abutment with Screw Brånemark System® NP	Brånemark System® MK III NP and NobelSpeedy® Groovy NP implants	Guided Sleeve NP		
Guided Template Abutment con Screw Brånemark System® RP	Brånemark System® MK III RP and NobelSpeedy® Groovy RP implants	Guided Sleeve RP		
Guided Template Abutment with Screw Brånemark System® WP	Brånemark System® MK III WP and NobelSpeedy® Groovy WP implants	Guided Sleeve 6.0/ WP		
Guided Template Abutment with Screw NobelActive® WP 5.5	NobelActive®/ NobelActive® TiUltra™ WP 5.5			
Handle for Guided Drill Guide	N/D	N/D	N/D	Guided Drill Guide

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Nome prodotto	Connessioni al manipo
Guided Start Drill/Counterb Bmk Syst NP	EN ISO 1797
Guided Start Drill/Counterbore NobRpl NP	
Guided Start Drill/Counterb Bmk Syst RP	
Guided Start Drill/Counterbore NobRpl RP	
Guided Start Drill/Counterb BmkSyst WP Ø 5	
Guided Start Drill/Counterbore NobRpl WP	
Guided Start Drill/Counterb BmkSyst WP Ø 6	
Guided Start Drill/Counterb NobRpl 6.0	
Guided Twist Drill 2x(10+)7-18 mm	
Guided Twist Drill 2.8x(10+)7-18 mm	
Guided Twist Drill 3x(10+)7-18 mm	
Guided Twist Drill 3.2x(10+)7-18 mm	
Guided Twist Drill 3.4x(10+)7-18 mm	
Guided Twist Drill 3.8x(10+)7-18 mm	
Guided Twist Drill 4.2x(10+)7-18 mm	
Guided Twist Drill 5x(10+)7-18 mm	
Guided Dense Bone Drill Taper NP 3.5x13 mm	
Guided Dense Bone Drill Taper NP 3.5x16 mm	
Guided Dense Bone Drill Taper RP 4.3x13 mm	
Guided Dense Bone Drill Taper RP 4.3x16 mm	
Guided Dense Bone Drill Taper WP 5x13 mm	
Guided Dense Bone Drill Taper WP 5x16 mm	
Guided Dense Bone Drill Taper 6.0 6x13 mm	
Guided Dense Bone Drill Taper 6.0 6x16 mm	
Guided Screw Tap NP 10-13 mm	
Guided Screw Tap RP Ø 3.75 7-13 mm	
Guided Screw Tap Tapered RP	
Guided Screw Tap WP Ø 5 7-13 mm	
Guided Screw Tap Tapered WP	
Guided Screw Tap 6.0/WP Ø 6 7-13 mm	
Guided Screw Tap Tapered 6.0	
Guided StartDrill/Countb former MkIII RP	
Guided Screw Tap RP Ø 4 7-13 mm	
Guided Twist Drill 3.2x(10+)7-13 mm	
Guided Twist Drill 3.4x(10+)7-13 mm	
Guided Twist Drill 3.8x(10+)7-13 mm	
Guided Twist Drill 4.2x(10+)7-13 mm	
Guided Twist Drill 5x(10+)7-13 mm	
Guided Twist Step Drill Ø 2.4/2.8 (10+) 7-13 mm	
Guided Twist Step Drill Ø 2.8/3.2 (10+) 7-13 mm	
Guided Twist Step Drill Ø 3.2/3.6 (10+) 7-13 mm	
Guided Twist Step Drill Ø 3.8/4.2 (10+) 7-13 mm	
Guided Twist Step Drill Ø 4.2/4.6 (10+) 7-13 mm	
Guided Twist Step Drill Ø 2.4/2.8 (10+) 7-18 mm	
Guided Twist Step Drill Ø 2.8/3.2 (10+) 7-18 mm	
Guided Twist Step Drill Ø 3.2/3.6 (10+) 7-18 mm	
Guided Twist Step Drill Ø 3.8/4.2 (10+) 7-18 mm	
Guided Screw Tap NobelActive® NP 8.5-10 mm	
Guided Screw Tap NobelActive® RP 4.3 8.5-10 mm	
Guided Screw Tap NobelActive® RP 5.0 8.5-10 mm	

Guided Twist Drill Tapered 2 x (10+) 8-16 mm
Guided Start Drill
Guided Counterbore NobelReplace® 3.5
Guided Counterbore NobelReplace® 4.3
Guided Counterbore NobelReplace® 5.0
Guided Counterbore NobelReplace® 6.0
Guided Screw Tap Tapered NP
Guided Pilot Twist Drill 1.5x(10+)10-15 mm
Guided Twist Step Drill 4.2/5.0 (10+) 7-13 mm
Guided Twist Step Drill 4.2/5.0 (10+) 7-18 mm
Guided Screw Tap NobelActive® WP 5.5 7-8.5 mm
Guided Screw Tap NobelActive® WP 5.5 10-15 mm
Guided Twist Step Drill Ø 4.2/4.6x(10+)7-18 mm
Guided Screw Tap NobelParallel™ CC 3.75 7-10 mm
Guided Screw Tap NobelParallel™ CC 3.75 11.5-13 mm
Guided Screw Tap NobelParallel™ CC 3.75 15-18 mm
Guided Screw Tap NobelParallel™ CC 4.3 7-10 mm
Guided Screw Tap NobelParallel™ CC 4.3 11.5-13 mm
Guided Screw Tap NobelParallel™ CC 4.3 15-18 mm
Guided Screw Tap NobelParallel™ CC 5.0 7-10 mm
Guided Screw Tap NobelParallel™ CC 5.0 11.5-13 mm
Guided Screw Tap NobelParallel™ CC 5.0 15-18 mm
Guided Screw Tap NobelParallel™ CC 5.5 7-10 mm
Guided Screw Tap NobelParallel™ CC 5.5 11.5-13 mm
Guided Screw Tap NobelParallel™ CC 5.5 11.5-15 mm
Guided Tissue Punch NP-sleeve (3.75)
Guided Tissue Punch WP-sleeve (5.5)
Guided Counterbore NobelParallel™ CC 3.75
Guided Counterbore NobelParallel™ CC 4.3
Guided Counterbore NobelParallel™ CC 5.0
Guided Counterbore NobelParallel™ CC 5.5
Guided Dense Bone Screw Tap NobelActive® NP 8.5-18 mm
Guided Dense Bone Screw Tap NobelActive® RP 4.3 8.5-18 mm
Guided Dense Bone Screw Tap NobelActive® RP 5.0 8.5-18 mm
Guided Dense Bone Screw Tap NobelActive® WP 5.5 7-8.5 mm
Guided Dense Bone Screw Tap NobelActive® WP 5.5 10-15 mm
Guided Drill Tapered NP
Guided Drill Tapered RP
Guided Drill Tapered WP
Guided Drill Tapered 6.0

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Compatibilità NobelProcera®/Procera® Abutment Titanium

Abutment in titanio - con Unigrip™

Connessione	Piattaforma	Vite clinica	Vite da laboratorio	Cacciavite	Torque
Connessione conica	3.0	Clinical Screw CC 3.0	Lab Screw Implant Level CC 3.0	Unigrip™	15 Ncm
	NP	Clinical Screw CC NP	Lab Screw CC NP		35 Ncm
	RP	Clinical Screw CC RP/WP	Lab Screw CC RP/WP		
	WP				
Trilobata	NP	Abutment Screw NobelReplace® NP	Lab Screw Implant Level NobRpl NP		
	RP	Abutment Screw NobelReplace® RP/WP/6.0	Lab Screw Implant Level NobRpl RP/WP/6.0		
	WP				
	6.0				
Esagonale esterna	NP	Abutment Screw Bmk Syst NP	Lab Screw CC NP		
	RP	Abutment Screw Bmk Syst RP	Lab Screw Ext Hex RP		
	WP	Abutment Screw Bmk Syst WP	Lab Screw Implant Level Bmk Syst WP		

Compatibilità NobelProcera®/Procera® Abutment Zirconia

Abutment in zirconio - con Unigrip™

Torque	Torque	Torque	Torque	Torque	Torque
Connessione conica	NP	Clinical Screw CC NP	Lab Screw CC NP	Unigrip™	35 Ncm
	RP	Clinical Screw CC RP/WP	Lab Screw CC RP/WP		
Trilobata	NP	Screw Ceramic Abutment NobRpl NP	Lab Screw Implant Level NobRpl NP		
	RP	Screw Ceramic Abutment NobRpl RP/WP/6.0	Lab Screw Implant Level NobRpl RP/WP/6.0		
	WP				
	6.0				
Esagonale esterna	NP	Screw Ceramic Abutment Bmk Syst NP	Lab Screw CC NP		
	RP	Screw Ceramic Abutment Bmk Syst RP	Lab Screw CC RP/WP		
	WP	Screw Ceramic Abutment Bmk Syst WP	Lab Screw Implant Level Bmk Syst WP		

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IFU3018. Ponte implantare in zirconia NobelProcera®

	Tipo di connessione											
	Esagonale esterna			Interna trilobata			CC			MUA		
NobelProcera Implant Bridge Zirconia	NP	RP	WP	NP	RP	WP	6.0	NP	RP	WP	NP/RP	WP
Adattatore clinico	N/D							Metal Adapter Implant Bridge Clinical CC			N/D	
Adattatore di laboratorio	N/D							Metal Adapter Implant Bridge Lab CC			N/D	
Vite clinica	Screw Ceramic Abutment Brånemark System®							Omnigrip™ Clinical Screw CC			Prosthetic Screw Multi-unit Abutment Omnigrip™ Mini	
Vite da laboratorio	Lab Screw Implant Level External Hex		N/D (*)				Omnigrip™ Lab Screw CC Lab Screw Implant Level Conical Connection			Lab Screw Multi-unit Abutment Omnigrip™ Mini		
Replica dell'impianto/MUA	Implant Replica Brånemark System®			Implant Replica NobelReplace®			Implant Replica Conical Connection			Multi-unit Abutment Plus Replica		Multi-unit Abutment Plus Replica
	Elos Accurate Analog for Printed Model External Hex			Elos Accurate Analog for Printed Model Tri-Channel			IOS Model Implant Replica CC			Elos Accurate Analog for Printed Model Multi-unit	Abutment Replica Multi-unit Bmk Syst Elos Accurate Analog for Printed Model Multi-unit	
Analoghi di protezione	Protection Analog Bmk Syst			Protection Analog NobelReplace®			Protection Analog/Drill Guide CC			Protection Analog Multi-unit	Protection Analog Multi-unit Protection Analog M-u Bmk Syst	
MUA	Multi-unit Abutment Bmk Syst 17° Multi-unit Abutment Bmk Syst 30° Multi-unit Abutment Bmk Syst Zygoma Abutment Brånemark System® Zygoma Zygoma 17° Brånemark System® Zygoma 17° 45° Multi-unit Abutment Ext Hex 60° Multi-unit Abutment Ext Hex			Multi-unit Abutment NobRpl 17° Multi-unit Abutment NobRpl 30° Multi-unit Abutment NobRpl			Multi-unit Abutment Plus Conical Connection 17° Multi-unit Abutment Plus Conical Connection 30° Multi-unit Abutment Plus Conical Connection Multi-unit Abutment Xeal™ Conical Connection 17° Multi-unit Abutment Xeal™ Conical Connection 30° Multi-unit Abutment Xeal™ Conical Connection			N/D		
Impianti	NobelSpeedy® Groovy NobelSpeedy® Shorty	Brånemark System® Mk III Brånemark System® Mk III TiUnite® NobelSpeedy® Groovy NobelSpeedy® Shorty Impianto Zygoma NobelZygoma™	NobelSpeedy® Groovy NobelSpeedy® Shorty	Replace Select™ TC			NobelActive® NobelActive® TiUltra™ NobelParallel™ CC NobelParallel™ CC TiUltra™ NobelReplace® CC NobelReplace® CC PMC NobelReplace® CC TiUltra™	NobelActive® NobelActive® TiUltra™ NobelParallel™ CC NobelParallel™ CC TiUltra™ NobelReplace® CC NobelReplace® CC PMC NobelReplace® CC TiUltra™	NobelActive® NobelActive® TiUltra™ NobelParallel™ CC NobelParallel™ CC TiUltra™	N/D		

IFU3020. NobelProcera® Zirconia N1 Base

Compatibilità NobelProcera® Zirconia NB N1™ Base – Vite protesica – Nobel Biocare N1™ Base Xeal TCC Tri

Tipo di protesi	Vite protesica	Cacciavite	Nobel Biocare N1™ Base Xeal™ TCC Tri
Nobel Procera® Zirconia Abutment/ Implant Crown NB N1™ Base NP	Prosthetic Screw NobelProcera® Zirconia Nobel Biocare N1™ Base NP	Omnigrip™ Mini	Nobel Biocare N1™ Base Xeal™ TCC Tri NP
Nobel Procera® Zirconia Abutment/ Implant Crown NB N1™ Base RP	Prosthetic Screw NobelProcera® Zirconia Nobel Biocare N1™ Base RP		Nobel Biocare N1™ Base Xeal™ TCC Tri RP

Compatibilità NobelProcera® Zirconia NB N1™ Base – Vite da laboratorio - Base Replica

Tipo di protesi	Vite da laboratorio	Base Replica Nobel Biocare N1™ Base	IOS Base Replica Nobel Biocare N1™ Base
Nobel Procera® Zirconia Abutment/ Implant Crown NB N1™ Base NP	Lab Screw NobelProcera® Zirconia Nobel Biocare N1™ Base NP	Base Replica Nobel Biocare N1™ Base Tri NP	IOS Base Replica Nobel Biocare N1™ Base Tri NP
Nobel Procera® Zirconia Abutment/ Implant Crown NB N1™ Base RP	Lab Screw NobelProcera® Zirconia Nobel Biocare N1™ Base RP	Base Replica Nobel Biocare N1™ Base Tri RP	IOS Base Replica Nobel Biocare N1™ Base Tri RP